HbA1c 2R CONTROL LOW

Cat. No.	Pack Name	Packaging (Content)
XSYS0098	HbA1c 2R CON L	4 × 0.5 ml

INTENDED USE

HbA1c 2R CON L is a set of four lyophilized controls based on human blood material (erythrocytes). HbA1c concentration in HbA1c 2R CON L is normal.

STORAGE

The controls both unopened and opened must be stored at 2–8 °C, protected from light and heat.

STABILITY

Unopened controls are stable until the date of expiry stated at the label.

After opening and reconstitution are controls stable for one week if stored under 2–8 °C and if contamination and evaporation are avoided after having opened the vials.

Proper storage and handling must be observed.

WARNING AND PRECAUTIONS

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent of the kit is not classified as dangerous.

Each individual blood donation used for production of HbA1c 2R CON L was found to be non-reactive when tested with approved methods for HBs-Ag, HIV-Ag+Ab and HCV-Ab HTLV-1-Ab. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the controls with the same precautions used for patient specimens.

PREPARATION

Before use perform reconstitution of control by adding 500 µl of hemolyzing solution (R3) into the vial, mix gently and allow to stand for 5 minutes. After reconstitution is control ready to use. Don't prepare hemolysate.

PROCEDURE

Please refer to the reagent package insert for instructions for use HbA1c 2R, Cat. No. XSYS0096.

ASSIGNED VALUE

Target values according to IFCC (mmol/mol):

	Lot No.	Expiry Date	Targer value (mmol/mol)	Range (mmol/mol)	1 SD (mmol/mol)
HbA1c Control Low	xxxxxxx	yyyy-mm	xx.x	xx.x–xx.x	x.x

Target values according to DCCT/NGSP (% NGSP):

	Lot No.	Expiry Date	Targer value (%)	Range (%)	1 SD (%)
HbA1c Control Low	xxxxxxx	yyyy-mm	x.xx	x.xx–x.xx	x.xx

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LITERATURA

- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329:977-86.
- Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
- Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Cem Lab Med 2002; 40:78- 89.
- Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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- ЛИТЕРАТУРА**
1. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329:977-86.
 2. Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
 3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78- 89.
 4. Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
 5. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
 6. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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- ЛІТЕРАТУРА**
1. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med. 1993; 329:977-86.
 2. Little RR, Rohlfing CI, Wiedmeyer HM, Myers GL et al. The national Glycohemoglobin Standardization Program: A five-Year Progress Report. Clin Chem 2001; 47:1985-92.
 3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC Reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002; 40:78- 89.
 4. Hoelzel W, Weykamp C et al. IFCC Reference system for measurement of HbA1c in Human Blood and the national standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study. Clin Chem 2004; 50:1:1666-74.
 5. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
 6. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).

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