

ADVISORY NOTICE

Update of target ranges for Bio-Rad Liquichek – batches: 87730 - 87760

LAURA XL, Catalogue No.: INS00065
Bio-Rad Liquichek Urinalysis Control, Catalogue No.: 435

Dear customers,

We would like to inform you, that according to our latest observation and customers feedback we updated target ranges of sediment particles (RBC and WBC) for Bio-Rad Liquichek Urinalysis Control – batches 87730 (exp: 02/2023), 87740 (exp: 07/2023), 87750 (exp: 12/2023) and 87760 (exp: 06/2024). The updated target values for sediment particles are summarized in the following table:

Bio-Rad Liquichek - batch 87730 (exp: 02/2023)		
Parameter	Level 1	Level 2
Red Blood Cells (RBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	24 - 189 N/hpf (132 - 1053 N/ μ l)
White Blood Cells (WBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	5 - 66 N/hpf (26 - 368 N/ μ l)
Bio-Rad Liquichek - batch 87740 (exp: 07/2023)		
Parameter	Level 1	Level 2
Red Blood Cells (RBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	20 - 199 N/hpf (105 - 1105 N/ μ l)
White Blood Cells (WBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	5 - 85 N/hpf (26 - 474 N/ μ l)
Bio-Rad Liquichek - batch 87750 (exp: 12/2023)		
Parameter	Level 1	Level 2
Red Blood Cells (RBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	24 - 199 N/hpf (132 - 1105 N/ μ l)
White Blood Cells (WBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	9 - 85 N/hpf (53 - 474 N/ μ l)
Bio-Rad Liquichek - batch 87760 (exp: 06/2024)		
Parameter	Level 1	Level 2
Red Blood Cells (RBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	14 - 199 N/hpf (79 - 1105 N/ μ l)
White Blood Cells (WBC)	0 - 3 N/hpf (0 - 16 N/ μ l)	5 - 85 N/hpf (26 - 474 N/ μ l)

In case of any questions, please contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.
 Technical Product Management
 Brno, 11th January 2023



Erba Lachema s.r.o.
 Karásek 2219/1d, 621 00 Brno
 DIČ: CZ26918846 -7-

