Uro-dip® 10 

REAGENT STRIPS FOR URINALYSIS / TESTSTREIFEN FÜR DIE URINUNTERSUCHUNG
TIRAS REACTIVAS PARA URIANÁLISIS / BANDETTES RÉACTIVES POUR LES ANALYSES D'URINE

	Specific gravity Spezifische Masse Peso específico Densité	Leucocytes Leukozyten Leucocitos Leucocytes Nitrites	Nitrites Nitrite Nitritos Nitrites	pH	Protein Eiweißstoffe Proteínas Protéine	Glucose Glukose Glucosa Glucose	Ketones Ketone Cetonas Cétones	Urobilinogen Urobilinogen Urobilínogeno Urobilinoène	Bilirubin Bilirubin Bilirubina Bilirubine	Blood Blut Sangre Sang
Uro-dip® 10e	URPH0021	100	■	■	■	■	■	■	■	■

(EN)

Reagent Strips for the rapid determination of Specific Gravity, Leucocytes, Nitrite, pH, Protein, Glucose, Ketones (Acetoacetic Acid), Urobilinogen, Bilirubin and Blood in urine. Diagnostic strips Uro-dip 10e are intended for objective evaluation by urine analysers Uro-dipcheck 240e and Uro-dipcheck 400e.

PRODUCT NAME: Uro-dip 10e

SUMMARY AND EXPLANATION: Uro-dip Reagent Strips are dip-and-read test strips for In Vitro Diagnostic Use only for testing above items in urine. Test result may provide information regarding the status of carbohydrate metabolism, kidney and liver function, acid-base balance, and urinary tract infection.

WARNING AND PRECAUTIONS: For in vitro diagnostic use only. For professional use only.

CHEMICAL PRINCIPLES OF PROCEDURE AND INGREDIENTS:

Specific gravity - The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in urine. Its result is colour change of acid-base indicator from blue-green colour in urine with low concentration of ions, through green and yellow-green in urine with increased concentration of ions, to amber yellow colour. Using this test it is possible to determine specific gravity of urine in the range of 1.000 up to 1.030. The first morning urine of healthy person should be in the range of 1.015 up to 1.025. **Ingredients:** poly(methylvinylether/maleic acid) 32 %; bromthymol blue 5.1 %

Leucocytes - The test is based on enzymatic reaction. The test pad contains an indoxyl ester, that is cleaved by granulocyte esterases. The released indoxyl reacts with a diazonium salt and pink or violet coloration is formed. The colour intensity is proportional to leucocytes amount in a sample of tested urine. **Ingredients:** indoxyl ester 0.43 %; diazonium salt 0.05 %

Nitrite - The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in urine. The colour test is based on the principle of the Griess test. Some degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10³ or more organisms in 1 ml of urine specimen, but coloration of pad is not quantitatively proportional to the amount of bacteria present in urine. Negative results do not exclude significant bacteriuria, as insufficient incubation may have occurred and some organisms causing urinary tract infections do not contain nitrate reductase to convert nitrate to nitrite. For these reasons the identification of known positive cases with the nitrite test is about 70%. We recommend to test the first morning urine specimens, when long bladder incubation has occurred. **Ingredients:** sulphanilamide 5.1 %; tetrahydrobenzo-[h]-quinoline 5.8 %

pH - The test is based on the double indicator principle and gives a range of colours from orange through yellow and green to blue and permits differentiation to within 0.5 pH unit in the range pH 5 to 9. **Ingredients:** methyl red 0.71 %; bromthymol blue 12.1 %

Protein - The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein. **Ingredients:** tetrabromphenolphthalein ester 0.21 %; tetrabromphenol blue 0.35 %

Glucose - The test is based on the specific glucose oxidase/peroxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars, it reacts with presence of D-glucose by greeny to dark-green coloration. **Ingredients:** glucose oxidase 1.3 %; peroxidase 1.3 %; tetramethylbenzidine 21 %

Ketones - The test is based on the principle of Legal's test and is more susceptible to acetoacetic acid than to acetone. Test does not react with β-hydroxybutyric acid. The colour scale is calibrated for acetoacetic acid. **Ingredients:** sodium nitroprusside 4.9 %

Urobilinogen - The test is based on the coupling of urobilinogen with stabilized reagent. The test is specific for urobilinogen and stercobilinogen and is not susceptible to the interfering factors known at Ehrlich's test. **Ingredients:** diazonium salt 2.3 %

Bilirubin - The test is based on the coupling of bilirubin with stabilized reagent. **Ingredients:** diazonium salt 0.75 %

Blood - The test is based on the peroxidase activity of haemoglobin which catalyzes the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales; for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin and may detect its presence from concentrations corresponding approx. to 5 Ery/μl. **Ingredients:** tetramethylbenzidine 1.5 %; cumene hydroperoxide 15.2 %

Compensation field – Pad, which is not impregnated with any reagents. The compensation field is used for suppressing of the dark colour of the urine sample, because the dark colour could have the effect on the evaluation of reagent pads.

STORAGE AND HANDLING:

Store in a cool, dry place at temperatures between (+2 to +30)°C. Do not store the strips in a refrigerator or freezer. Store away from moisture and light. When stored in a original container, the product is stable up to the expiry date printed on the label.

Replace the bottle cap immediately and tightly after removing test strips, and keep the vial tightly closed between tests.

Do not touch test areas of urine reagent strips. Do not open container until ready to use.

Discoloration or darkening of the test pads may indicate deterioration. If this is evident, or if test results are questionable or inconsistent with expected finding, confirm that the product is within its expiration date and is reacting properly using known negative and positive control materials.

SPECIMEN COLLECTION AND PREPARATION:

Collect urine in a clean, dry container that allows complete immersion of all the fields on the test strip. Do not add preservatives. Test the specimen as soon as possible, with the sample well mixed but not centrifuged. The use of fresh morning urine is recommended for optimal nitrite tests, as well as for the valid determination of bilirubin and urobilinogen, since these compounds are unstable when exposed to light. If immediate testing is not possible, the sample should be stored in the refrigerator, but not frozen, and then brought to room temperature before used in the test. Unpreserved urine at room temperature may undergo pH changes due to microbial proliferation, which may interfere with protein determination. If cleanly voided specimens are not collected from females, positive results for leucocytes may be found due to contamination from outside the urinary tract. Skin cleansers containing chlorhexidine may affect protein test results if specimen contamination occurs.

OBJECTIVE TEST PROCEDURE:

The procedure must be followed exactly to achieve reliable results. Do not divide test strips!

- 1) Remove only as many test strips as are required, and reseal the tube immediately after use.
- 2) Do not touch test pads of the strips.
- 3) Completely immerse all reagents pads in specimen (no longer than 1-2 sec.)
- 4) Run edge of the strip against rim of urine container to remove excess urine. Hold the strip in horizontal position.
- 5) Evaluate the result using the reader Uro-dipcheck® 240e or Uro-dipcheck® 400e, follow enclosed instructions for that instruments.

Notes: For visual evaluation compare the tests pads to the corresponding colour scale on the label after approx. 60 sec. and for leucocytes after 120 sec. According to the different spectral sensitivity of the human eye and the optical system it is not possible to guarantee the exact correspondence between visual reading and reading by analysers.

QUALITY CONTROL:

For best results, performance of reagent strips should be confirmed by testing known negative and positive specimen or controls (e.g. URINORM) whenever a new bottle is first opened. Each laboratory should establish its own goals for adequate standards of performance. Each lab worker should ensure that it complies with government and local requirements.

LIMITATIONS OF PROCEDURE:

Specific gravity - The reaction is not affected by pH values of urine over 6.5 shift colour response towards lower values of specific gravity.

Leucocytes - In case when the urine sample is more markedly cloudy (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of colour reaction is increased by alkaline pH and higher urine density.

Nitrites - Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy and vitamin C for 3 days prior to the test. Sensitivity of this test declines with high specific gravity of urine. Increased diuresis can cause false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of urine. Test can be applied only at fresh urine. Inaccurate results may occur at stale urines, in which nitrite can be formed by contamination of the specimen.

Protein - In strongly alkaline urines (pH >8) from patients on medication with quinine or quinine containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quaternary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false-negative results may occur. Do not take the colour of the dry pad into consideration.

Glucose - The reaction is independent on pH and presence of ketone bodies.

Ketones - Drugs and diagnostics on the basis of phenolphthalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.

Urobilinogen - The reaction is not affected by pH of urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.

Bilirubin - The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 μmol/l) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine). The reaction is not affected by pH of urine.

Blood - Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmacological origin.

All diagnostics pads do not interfere with the common concentration of the ascorbic acid.

EXPECTED VALUES:

Specific gravity - The normal SG of urine ranges from 1,015 to 1,025. / **Leucocytes** - < 10 Leu/μl / **Nitrites** - Normally no nitrites are detectable in urine.

pH - Normal values of pH lies between 5.5 to 7. / **Protein** - Normal urine specimens ordinarily contain some protein (< 15 mg/dl).

Glucose - Abnormal concentration of glucose in urine is higher than 25 mg/dl. / **Ketones** - Normal concentration of ketone bodies is lower than 2 mg/dl.

Urobilinogen - Normal urobilinogen concentration in urine is lower than 1 mg/dl. / **Bilirubin** - Normal bilirubin concentration in urine is lower than 0,2 mg/dl.

Blood - < 5 Ery/μl.

TEST PAD AND SENSITIVITY (SPECIFICITY):

Leucocytes - 10 Leu/μl - intact and lysed WBCs / **Nitrites** - 11 mmol/l (0,05 mg/dl) - nitrites / **Protein** - 0,15 g/l (15 mg/dl) - albumin

Glucose - 0,9 mmol/l (16 mg/dl) - D-glucose / **Ketones** - 0,1 – 0,2 mmol/l (1,0 – 2,0 mg/dl) - acetoacetic acid

Urobilinogen - 6,0 mol/μl (0,35 mg/dl) - urobilinogen, stercobilinogen / **Bilirubin** - 4,3 – 5,2 μmol/l (0,25 – 0,30 mg/dl) - bilirubin

Blood - 5 Ery/μl - haemoglobin

PLEASE NOTE: Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after discontinuing a drug. The sensitivity depends upon the variability of urines. The semiquantitative analysis is not sufficient for the completing of diagnoses.

QUALITY CONTROL:

URINORM control urines, Cat. No. REG00053 are designed for verification and confirmation of precision and accuracy of Uro-dip® 10e diagnostic strips as well as Uro-dipcheck 240e and Uro-dipcheck 400e. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in Uro-dipcheck 240e and Uro-dipcheck 400e user manuals.

STORAGE: Keep Diagnostic Test Strips in tightly closed original tubes in a dry and dark place at (+2 to +30)°C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack.

WASTE DISPOSAL: Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

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