BILIRUBIN DIRECT DCA

Cat. No.	Pack Name	Packaging (Content)
XSYS0086	BIL D DCA 330	R1: 6 x 44 ml, R2: 6 x 11 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Bilirubin Direct in human serum and plasma.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of haemoglobin. Bilirubin formed in the reticulo endothelial system is transported bound by albumin to the liver. This bilirubin is water insoluble and is known as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid to form direct bilirubin. Conjugated bilirubin is excreted via the biliary system into the intestine. Here it is metabolised by bacteria to urobilinogen and stercobilinogen.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

Total Bilirubin is elevated in obstructive conditions of the bile duct, hepatitis, cirrhosis, in haemolytic disorders and several inherited enzyme deficiencies.

Indirect Bilirubin is elevated by pre-hepatic causes such as haemolytic disorders or liver diseases resulting in impaired entry transport or conjugation within the liver. Monitoring of indirect bilirubin in neonates is of special importance as it is the indirect (or free) bilirubin bound to albumin that is able to cross the blood brain barrier more easily increasing the danger of cerebral damage.

PRINCIPLE

Conjugated (direct) bilirubin reacts with diazotized 2,4-dichloroaniline in acidic solution to produce an intensely coloured red diazo compound (520–560 nm). The intensity of color of this dye in solution is proportional to the concentration of direct bilirubin.

REAGENT COMPOSITION

R1	
EDTA	0.1 mmol/l
NaCl	0.26 mol/l
R2	
EDTA	0.1 mmol/l
Diazotized 2,4-dichloroaniline	0.1 mmol/l
HCl	0.18 mol/l

REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

After first opening, reagents are stable for 60 days at 2–8 °C if stored at appropriate conditions, closed carefully, without any contamination and away from light sources.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability: 2 days at 15–25 °C
7 days at 2–8 °C
3 months at -20 °C
FREEZE ONLY ONCE!

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to Jendrassik Grof original formulation.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

mg/dl x 16.95 = μmol/l

EXPECTED VALUES ³

Adults and infants : 0–0.2 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.064 mg/dl
Linearity: 14.75 mg/dl
Measuring range: 0.064–14.75 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	0.385	0.019	4.99
Sample 2	0.474	0.027	5.60

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	0.729	0.032	4.41
Sample 2	2.343	0.075	3.20

COMPARISON

A comparison between XL-Systems BIL D DCA (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.064 x – 0.011 mg/dl

r = 0.997

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, triglycerides up to 500 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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