Erba H7 ret CON

Haematology Control

Cat. No.	Product Name	Level	Packaging (Content)
HEM00059	Erba H7 ret CON	Low, Normal, High	3 × 3 mL, Instruction for use
HEM00060	Erba H7 ret CON L	Low	1 × 3 mL, Instruction for use
HEM00061	Erba H7 ret CON N	Normal	1 × 3 mL, Instruction for use
HEM00062	Erba H7 ret CON H	High	1 × 3 mL, Instruction for use

INTENDED USE

Erba H7 ret CON is an *in vitro* diagnostic reagent used for control of the accuracy and precision of RET measurements on automated haematology analysers H7100. For professional use in clinical laboratory only.

PRINCIPLE

It is an established laboratory practice to use a stable control to monitor the performance of diagnostic tests. This control is composed of stable materials that provide a means of monitoring the performance of haematology blood cell counters. It is sampled in the same manner as a patient specimen.

DESCRIPTION AND COMPOSITION

Erba H7 ret CON is an *in vitro* diagnostic material composed of human and porcine cells suspended in a plasma-like fluid with preservatives.

MEASURED PARAMETERS

Parameter	Abbreviation	Default Unit
Number of Reticulocyte	RET	10 ¹² /L
Percentage of Reticulocyte	RET%	%

CALCULATED PARAMETERS

Parameter	Abbreviation	Default Unit
Low Fluorescent Ratio	LFR	%
Middle Fluorescent Ratio	MFR	%
High Fluorescent Ratio	HFR	%
Immature Reticulocyte Fraction	IFR	%
Reticulocyte Hemoglobin Expression	RHE	pg

PREPARATION

The haematology controls are ready to use.

STABILITY AND STORAGE

Store Erba H7 ret CON upright at 2–8 °C when not in use. Protect tubes from overheating and freezing. Unopened tubes are stable through the expiration date. Opened tubes are stable for 15 days at 2–8 °C if they are handled properly.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser: H7100, Cat. No. INS00087

ASSAY PROCEDURE

- Remove tubes from the refrigerator and allow them to warm to room temperature (15–25 °C) for 15 minutes before mixing.
- To mix, hold a tube horizontally between the palms of the hands. Do not premix on a mechanical mixer.
 - Roll the tube back and forth for 20–30 seconds; occasionally invert the tube. Mix vigorously, but do not shake.
 - Continue to mix in this manner until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.
 - Gently invert the tube 8–10 times immediately before sampling.
- Analyze the sample as instructed in the Quality Control section of the User Manual for your instrument.
- After sampling:
 - If tube has been open for sampling, clean residual material from the cap and tube rim with a lint-free tissue. Replace the cap tightly.
 - Return tubes to refrigerator within 30 minutes of use.

Note: After mixing, product should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. Do not use the product if deterioration is suspected.

EVALUATION

The H7100 analyser performs the evaluation automatically.

EXPECTED RESULTS

You can find the assay values sheet on the website <https://www.erbalachema.com/en-hematology-files/>. Search the assay values according to the name of the lot. Verify that the lot number on the tube matches the lot number on the table of assay values. Assay values are determined on well-maintained, properly calibrated instruments using the instrument manufacturer's recommended reagents. Reagent differences, maintenance, operating technique, and calibration may contribute to inter-laboratory variation.

Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.

Assay values of a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument is in good working order and quality

control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.

For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the assay sheet if the control is suitable for the method.

Erba Lachema recommends running the QC program on a daily basis with low, normal and high level controls. Perform the calibration, if the results of control bloods are out of the limitation.

TRACEABILITY

The controls are traceable to the internationally agreed reference method.

- WBC:** A series of 1:500 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
- RBC:** A series of 1:50 000 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
- HGB:** Haemoglobin value is determined by spectrophotometric procedure according to CLSI Standard H15-A3 and is traceable to ICSH/WHO International Haemoglobin-cyanide Standard.
- HCT:** Packed cell volume (PCV) is measured by the microhaematocrit procedure according to CLSI Standard H07-A3. No correction is made for trapped plasma.
- MCV:** MCV is the calibrated parameter instead of the HCT. The MCV is calculated from the HCT and RBC using the formula: $MCV = HCT \times 10/RBC$.
- PLT:** A series of 1:126 dilutions are made using calibrated glassware in 1 % ammonium oxalate. Platelets are counted using a haemocytometer and phase contrast microscopy.

PERFORMANCE CHARACTERISTICS

All performance characteristics of the H7100 haematology systems are listed in the H7100 manual in the specification chapter.

LIMITATIONS

The performance of this product is assured only if it is properly stored and used as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

POTENTIAL BIOHAZARDOUS MATERIAL

Each human donor/unit used in the preparation of this product has been tested and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of tubes follow precautions for patient specimens as specified in biosafety procedures.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Erba H7 ret CON should be disposed as potentially infectious material. Please refer to local legal requirements.

REFERENCES

- Clinical Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood: Approved Standard-Third Edition. CLSI document H15-A3.
- Clinical Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhaematocrit Method: Approved Standard-Third Edition. CLSI document H07-A3.

USED SYMBOLS

 Catalogue Number

 Lot Number

 Expiry Date

 Contains human blood or plasma derivatives

 *In vitro* diagnostic medical device

 Manufacturer

 Storage Temperature

 Contains biological material of animal origin

 See Instructions for Use

 Content

 Biological Risks

 Contains biological material of human origin

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REFERENCE

- Clinical Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood: Approved Standard-Third Edition. CLSI document H15-A3.
- Clinical Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhematocrit Method: Approved Standard-Third Edition. CLSI document H07-A3.