Erba H7100 Clean Haematology Analyser Cleaner

Cat. No.	Pack Name	Packaging (Content)
HEM00063	Erba H7100 Clean 10 × 4 mL	10 × 4 mL, Instruction for use
HEM00064	Erba H7100 Clean 1 × 4 mL	1 × 4 mL, Instruction for use

INTENDED USE

Erba H7100 Clean is an *in vitro* diagnostic solution used for cleaning, rinsing and washing of automated haematology analysers H7100. For professional use in clinical laboratory only.

PRINCIPLE

The product is a strong alkaline solution that can effectively clean out residual blood components (deposited proteins, etc.) and other particles that affects blood cell count.

DESCRIPTION AND COMPOSITION

Sodium Hypochlorite; Sodium Hydroxide; Surfactant.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations.

Before using, leave the reagent at room temperature for at least 24 hours.

For reagent connection to the analyser, refer to the H7100 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser H7100, Cat. No. INS00087

SPECIFICATION

The product is a clear liquid without sediment, particles or floccules. Background result: WBC $\leq 0.1 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB $\leq 1 \text{ g/L}$, PLT $\leq 5 \times 10^9/L$, HCT $\leq 0.5 \%$.

STABILITY AND STORAGE

The product is stable till the expiry date mentioned on the label when stored at (+2 to +30) °C. Discard opened container after 60 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: JGQQ-6GE0-A35E-Q3D9

Danger

Reagent contains: sodium hypochlorite

Hazard statements:

H315 Causes skin irritation.

H318 Causes serious eye damage.

H412 Harmful to aquatic life with long lasting effects.

Precaution statements:

P280 Wear eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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