

## EC cartridge VET

Cat. No.	Pack Name	Number of Tests
VET00002	EC cartridge VET	250

**VETERINARY USE ONLY**  
**VETERINARY TECHNICAL DEVICE**



### INTENDED USE

Automatic analyser intended for *in vitro* quantitative determination of specific electrolytes in whole blood, serum, plasma of animals (see table below for full list). The instrument is able to determine the concentration of sodium (Na<sup>+</sup>), potassium (K<sup>+</sup>), ionized calcium (iCa<sup>2+</sup>) and chloride (Cl<sup>-</sup>). The results for tested electrolytes can be used for screening, monitoring and diagnosis of diseases associated with ionic imbalance as kidney disease, disorders of acid-base balance, nutritional disorders and digestive tract function. For professional use in clinical and veterinary laboratory only.

### PRINCIPLE

The electrolyte sensors are potentiometric sensors which measure the potential difference between an ion-selective electrode (ISE) and a reference electrode immersed in a sample liquid. The electrochemical potential between the sample side and the solid contact of the ion-selective membrane is proportional to the concentration of the corresponding ion in the sample. These potentials are recorded, electronically amplified and calculated<sup>1</sup>.

### DESCRIPTION AND COMPOSITION

#### 1) EC cartridge VET

##### Calibration solution 1:

Na<sup>+</sup> 150 mmol/L, K<sup>+</sup> 5 mmol/L, iCa<sup>2+</sup> 1.25 mmol/L, Cl<sup>-</sup> 103 mmol/L, buffer, preservative, wetting agent

##### Calibration solution 2:

Na<sup>+</sup> 105 mmol/L, K<sup>+</sup> 2 mmol/L, iCa<sup>2+</sup> 2.34 mmol/L, Cl<sup>-</sup> 71.1 mmol/L, buffer, preservative, wetting agent

##### Reference solution:

K<sup>+</sup> 1.2 mol/L, buffer, preservative, wetting agent

#### 2) EC Biosensor: Cat. No. REG00064 (for ordering extra pieces)

### PREPARATION

Insert EC Biosensor into the cartridge according to the EC 90 VET user manual.

### MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

EC 90 VET analyser, Cat. No. VET00001

Erba Clean I. Cat. No. EHL00042

Control material: ERBA NORM 4 × 5, Cat. No. BLT00080

ERBA NORM 10 × 5, Cat. No. XSYS0123

ERBA PATH 4 × 5, Cat. No. BLT00081

ERBA PATH 10 × 5, Cat. No. XSYS0124

### STABILITY AND STORAGE

Store at 2–30 °C. Note the expiration date on the labels.

On board stability: 90 days after opening. On-board stability of EC Biosensor is at least 45 days.

### SPECIMEN COLLECTION AND HANDLING

Whole blood samples may be directly measured or further processed to create plasma and serum samples. Take care to prevent haemolysis in blood samples during handling and processing of specimens. Use only sample containers spray-coated with sodium heparin. Do not use EDTA or citrate as anticoagulants. This may result in incorrect measurement results.

### QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in the EC 90 VET User manual.

### METROLOGICAL TRACEABILITY

Traceability is ensured by testing laboratory standards according to standard EN ISO 17511:2021 against primary and secondary certified reference materials.

### ASSAY PROCEDURE AND CALCULATION

- In the analyser software, in the maintenance screen, click on CHANGE EC CARTRIDGE button. Follow instructions on the screen.
- Unpack new cartridge, remove the cap and seal cap, insert EC Biosensor and place it in the analyser.
- In the analyser software click the OK button. The software will detect new cartridge and perform priming automatically.
- Calibration will be performed automatically.

The EC 90 analyser performs the evaluation automatically according to Nernst equation<sup>1</sup>:

$$E_i = E_i^0 + \frac{R \cdot T}{Z \cdot F} \ln a_i$$

Where:

$E_i$  Ion activity dependent potential

$E_i^0$  Standard potential

$R$  universal gas constant

$T$  absolute Temperature

$Z$  valency of ion  $i$

$F$  Faraday constant

$a_i$  activity of ion  $i$

### EXPECTED VALUES

The EC 90 VET has 13 different animals predefined – see the table below. For some of these animals, expected values are set according to the The Merck Veterinary Manual, 11th Edition<sup>2</sup>. Other animals (bird, rat, mouse, reptile, camel, other) have a *default expected values range*. However, all ranges can be adjusted – please see the User manual of the instrument:

	Na [mmol/L]	K [mmol/L]	iCa [mmol/L]	Cl [mmol/L]
cat	143–156	3.5–5.2	1.07–1.47	108–128
dog	142–152	3.9–5.1	1.18–1.37	110–124
rabbit	145–155	2.7–5	1.35–1.75	96–122
horse	137–148	2.8–5.1	1.25–1.65	99–109
cow	143–151	4.1–5.3	1.0–1.4	97–111
pig	135–150	4.4–5.1	0.89–1.45	94–106
sheep	135–145	3.5–5.1	1.44–1.6	98–107
goat	135–145	3.5–5.1	1.12–1.48	98–107
Default	135–145	3.5–5.1	1.15–1.3	98–107

### ANALYTICAL PERFORMANCE

	Specified for	Measuring range [mmol/L]	Display range [mmol/L]	Precision within run	Precision between runs	Resolution [mmol/L]
Na <sup>+</sup>	Whole blood Serum Plasma Aqueous	100–200	20–200	≤ 1 % CV	≤ 2 % CV	0.1
K <sup>+</sup>	Whole blood Serum Plasma Aqueous	1–15	0.2–4	≤ 1 % CV	≤ 2.5 % CV	0.01
iCa <sup>2+</sup>	Whole blood Serum Plasma Aqueous	0.1–2.5	0.1–6	≤ 0.02 SD in mmol/L	≤ 0.03 SD in mmol/L	0.01
Cl <sup>-</sup>	Whole blood Serum Plasma Aqueous	50–170	25–200	≤ 1 % CV	≤ 2.5 % CV	0.1

### INTERFERENCES

The following substances were found to significantly interfere with respective measurement parameters of the instrument.

Substance	Test concentration	Significant interference with				Test matrix
		c Na <sup>+</sup>	c K <sup>+</sup>	c iCa <sup>2+</sup>	c Cl <sup>-</sup>	
Magnesium (Mg <sup>2+</sup> )	> 2.2 mmol/L	–	–	✓	–	Plasma
Bromide (Br <sup>-</sup> )	> 5.3 mmol/L	–	–	–	✓	Plasma
Iodide (I <sup>-</sup> )	> 3.0 mmol/L	–	–	–	✓	Plasma
Thiocyanate (SCN <sup>-</sup> )	> 3.9 mmol/L	–	–	–	✓	Plasma
pH	6.9–7.9 pH units	–	–	✓	–	Plasma
		–	–	✓	–	Whole Blood
Citrate	> 0.4 mmol/L	–	–	✓	–	Plasma
EDTA	> 0.5 mmol/L	–	–	✓	–	Plasma

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

#### Calibration solution 1 and 2:

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

#### Reference solution:

The reagent is not classified as dangerous.


### WASTE MANAGEMENT


Users should be exceedingly careful when working with this biological material and treat the waste as potential infectious material. Please refer to local legal requirements.


### REFERENCES


- Orna, Mary Virginia; Stock, John (1989). Electrochemistry, past and present. Columbus, OH: American Chemical Society. ISBN 978-0-8412-1572-6. OCL 19124885.
- AIELLO, Susan E.; MOSES, Michael A.; ALLEN, Dana G.; CONSTABLE, Peter D.; DART, Andrew et al. (ed.). The Merck veterinary manual. Eleventh edition. Kenilworth, NJ: Merck & Co., 2016. ISBN 978-0-911-91061-2.


## USED SYMBOLS


 Catalogue number


 Lot number


 Expiry date


 Consult instructions for use

 Biological risks

 Manufacturer

 Temperature limit

 Content

 Contains sufficient for <n> tests