# EN-COCCUStest

Kat. č. MLT00004

## Pro mikrobiologii

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(FACKLAM, R. R. and COLLINS, M. D. 1989:

Identification of *Enterococcus* species isolated from human infections by a conventional test scheme.

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Cat. No.: MLT00004

## For microbiology

The kit EN-COCCUStest is designed for the routine identification of clinically significant representatives of *Enterococcus* family. The kit enables the identification of 36 strains by means of eight biochemical tests: arginine, sorbose, arabinose, sorbitol, melibiose, raffinose and mellezitose. For the detection of the characteristic feature of *Enterococcus* family – the pyrrolidonylarylamidase activity can be used PYRAtest, supplied in form of detection strips.

### The kit EN-COCCUStest contains:

- 3 microtitration plates (for identification of 12 strains each) with desiccant
- Instructions for use including the differentiation table
- Colour scale for EN-COCCUStest kit
- 3 polyethylene bags for incubation
- Storage bag (for storage of an open plate), 1 pc
- 36 record sheets
- Lid

### Storage, Shelf life:

EN-COCCUStest should be stored at the temperature of (+2 to +25)°C. Expiration is indicated on each package.

## Recommended working procedure for EN-COCCUStest

### Required material for performing the identification

(not included in the kit):

- Paraffin oil sterilized (Cat. No. MLT00042 – 750 determinations)
- Petri dishes with the blood agar
- Test tubes (100x15) mm with 1.2 ml of a sterile physiological saline
- Instrument DENSILAMETER II, Cat. No.: INS00062
- Automatic micropipette 0.1 ml, sterile tips
- Thermostat (35–37)°C
- Usual microbiological laboratory equipment (loops, markers, burner)

### For results evaluation (not included in the kit):

- Code Book for EN-COCCUStest - located at [www.erbalachema.com](http://www.erbalachema.com)
- The ErbaExpert Identification Program

### For detection of pyrrolidonylarylamidase activity:

- PYRAtest (Cat. No. MLT00040 – 50 determinations)
- Reagent for PYRAtest (Cat. No. MLT00023)

### Caution:

- For professional use only

### Follow the principles for working with infectious material!

### Isolation of cultures:

- Perform the isolation of cultures by usual technique on blood agar.
- Evaluate the morphology of the pure culture, the haemolytic activity on blood agar, perform test for the detection of the pyrrolidonylarylamidase activity, resp. other tests confirming the pertinence to *Enterococcus* family (usually serological group D or occasionally serological group Q or without any group antigene).

### Inoculum preparation:

- Prepare a suspension in the physiological saline from a pure, 24 hours old culture. Homogenize the suspension well.
- The suspension turbidity has to match a McFarland of 2.

#### Note:

- Usage of sampling pads is recommended to obtain correct suspension density

### Culture purity control:

- If required confirm purity of the suspension by streaking-out a sample from the inoculated physiological saline on blood agar and incubate for 24 hours. This culture can be used also for performance of complementary tests.

### Preparation of the EN-COCCUStest plate:

- Open the aluminium sachet close to the weld and take out the plate.
- Cut off required number of strips from the plate.
- Remove the adhesive tape from individual strips and insert them into prepared frame. In case you work with MIKRO-LA-TEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.
- Record number of the strains or isolates to be examined on the appropriate strips.
- Put the rest of the plate with desiccant to the storage bag and store it in a refrigerator for further use; protect it against humidity. We recommend to use the plate within 4 weeks after opening.

## Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

## Inoculation:

- Homogenize the bacterial suspension thoroughly before use.
- Inoculate 0,1 ml of the suspension into all microwells.
- After the inoculation, add 2 drops of the sterile paraffin oil into the wells of the column H (test arginine).

## Note:

- Inoculate control strains simultaneously with any doubtful strains or using a new batch of EN-COCCUStest to control colour expression of positive and negative reactions. To check negative reactions, drop the sterile physiological saline into one row.

## Note:

Clean the inside of the lid and the frame by ethanol just before the use.

## Incubation:

- Insert the frame with inoculated rows into the polyethylene bag for incubation.
- Fold the open end of the bag under the plate to prevent evaporation.
- Incubate at (+35 to +37)°C for 24 hours.

## Results reading:

- Read all tests and record the results in the record sheets.
- Use the table "Interpretation of reactions", the colour scale or colour reactions of the control strains to evaluate the colour reactions properly.

## Identification:

- Use Differentiation table or Index.
- Evaluate the culture as a complex considering the the origin of the isolate, the morphological features, the formation of the yellow pigment and the motion.
- If you have failed to identify the culture, repeat the procedure as above.

## Disposal of the material used:

- Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.
- Put paper packaging waste to recycling.

## The most frequent causes of identification failure:

- A mixed or contaminated culture.
- Inoculum of low density or small volume.
- Inoculum has contaminated adjacent strips.
- The test arginine has not been overlayed by paraffin oil.
- Failure to follow the recommended working procedure.
- It might be an atypical strain or representative of species that has not been included in the Differentiation table and Index.

## Identification by means of Index:

- The Index is intended for quick and standard interpretation of results. It includes typical results of most often isolated enterococci from a clinical material (the species included in the Index are marked in the Differentiation table with bold letters).
- Both the results and interpretations are arranged in the Index according to the profile values, i.e. the numerical code. The record sheet is designed to make the profile calculation.

## Procedure when calculating the profile:

- The results in the record sheet are subdivided into three groups.
- In each group, the tests are marked with numbers 1, 2 and 4.
- When calculating the profile, sum up the values of numbers with tests giving positive reactions in each group of tests.
- Record the sum into the respective field of the row "Profile" in the record sheet.
- Profile is represented by three digits number.

## Example of the Profile calculation:

H ARG 1	G SOE 2	F ARA 4	E MAN 1	D SOR 2	C MLB 4	B RAF 1	A MLZ 2	
+	—	—	+	+	—	—	—	— EN-COCCUStest
1	(0)	(0)	1	2	(0)	(0)	(0)	— numerical test values
1	1	1	1	2	1	1	1	— reaction results
1	(0)	(0)	1	2	(0)	(0)	(0)	— sum of values
1	1	1	1	3	1	1	1	— profile

The identification found in the Index, corresponding to the profile 130: *Enterococcus faecalis*.

**Quality control of tests:**

The quality control of the kits is performed systematically at various stages of their production. The batches are checked by tests on standard bacterial cultures. For those who wish to perform their own quality control tests, cultures mentioned in the table **Control strains** are recommended. **These strains are used to check the functionality of the kit, not to check the accuracy or success of the identification!**

- *Enterococcus raffinosus* CCM 4216 (ATCC 49427)
- *Enterococcus durans* CCM 5612 (ATCC 19432)

These strains are supplied by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz). The strains are delivered in a lyophilized form or on gelatinized disks.

**Control strains**

CCM No.	Column/Abbreviation of the test							
	H ARG	G SOE	F ARA	E MAN	D SOR	C MLB	B RAF	A MLZ
4216	–	+	+	+	+	+	+	+
5612	+	–	–	–	–	–	–	–

**Interpretation of reactions**

Column	Test	Abbreviation of the test	Reaction	
			positive	negative
H	Arginine	ARG	red-to-violet, red	yellow, pale orange
G	Sorbose	SOE	yellow, pale orange	red, orange-to-red
F	Arabinose	ARA	yellow, pale orange	red, orange-to-red
E	Mannitol	MAN	yellow, pale orange	red, orange-to-red
D	Sorbitol	SOR	yellow, pale orange	red, orange-to-red
C	Melibiose	MLB	yellow, pale orange	red, orange-to-red
B	Raffinose	RAF	yellow, pale orange	red, orange-to-red
A	Mellezitose	MLZ	yellow, pale orange	red, orange-to-red

**Health protection:**

Components of the kit are not classified as dangerous.

**Differentiation table**

V. 1.0

EN-COCCTest								Identification	Additional characters								
H ARG	G SOE	F ARA	E MAN	D SOR	C MLB	B RAF	A MLZ		PIG	MOT	SUC	RIB	LAC	TEL	SSD	PYR	HEM
+	-	-	+	+	-	-	(+)	<i>Enterococcus faecalis</i>	-	-	+	+	+	+	+	+	β 6%
+	-	+	+	(-)	(+)	(-)	-	<i>Enterococcus faecium</i>	-	-	(+)	+	+	-	+	+	β 0%
+	-	+	+	d	+	(+)	-	<i>Enterococcus mundtii</i>	+	-	+	+	+	-	+	+	
-	+	+	+	+	+	+	+	<i>Enterococcus raffinosus</i>	-	-	+	+	+	-	+	+	
-	+	+	+	+	-	-	+	<i>Enterococcus avium</i>	-	-	+	+	+	-	d	+	
d	-	+	+	d	+	+	-	<i>Enterococcus casseliflavus</i>	+	+	+	+	+	-	+	+	α
+	-	+	+	-	+	+	-	<i>Enterococcus gallinarum</i>	-	+	+	+	+	-	+	+	
+	-	-	-	-	d	-	-	<i>Enterococcus durans</i>	-	-	-	+	+	-	+	+	β 70%
+	-	-	-	-	(+)	d	-	<i>Enterococcus hirae</i>	-	-	(+)	+	+	-	d	+	β 0%
+	-	-	-	-	-	-	d	<i>E. faecalis asaccharolyt var.</i>	-	-	-	+	+	+	+	+	
-	+	-	+	+	(+)	+	-	<i>Enterococcus malodoratus</i>	-	-	+	+	+	-	d	+	
-	+	-	+	+	-	-	-	<i>Enterococcus pseudoavium</i>	-	-	+	+	+	-	d	+	
-	-	-	d	d	d	+	d	<i>Enterococcus cecorum</i>	-	-	+	+	+	n	-	-	
-	-	(+)	+	(+)	d	+	(-)	<i>Enterococcus columbae</i>	-	-	n	+	+	n	-	-	
+	-	-	-	-	-	+	-	<i>Enterococcus dispar</i>	-	-	+	+	+	-	-	+	
-	(-)	-	+	+	+	+	+	<i>Enterococcus saccharolyticus</i>	-	-	+	+	+	n	-	-	
+	-	-	+	+	(-)	-	+	<i>Enterococcus solitarius</i>	-	-	+	n	-	n	+	+	
-	-	-	-	-	d	+	+	<i>Enterococcus sulfureus</i>	+	-	+	+	+	-	-	+	

**Legend:**

+	= positive reaction	-	= negative reaction	d	= variable reaction
(+)	= prevailing positive reaction	(-)	= prevailing negative reaction	n	= unknown reaction
PIG	= formation of yellow pigment	RIB	= ribose	SSD	= serological group
MOT	= motion	LAC	= lactose	PYR	= pyrrolidonylarylamidase
SUC	= sacchare	TEL	= tellurite D	HEM	= haemolysis

*E. durans*, *E. hirae*, *E. faecalis asaccharolytic var.* belong to the group III *Enterococcus*

(FACKLAM, R. R., and COLLINS, M. D. 1989:

Identification of *Enterococcus species* isolated from human infections by a conventional test scheme. J. Clin. Microbiol. 27, 731-734).

The strains that could not be identified as belonging to any of the three species mentioned, should be identified as members of the *Enterococcus* Group III.

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