ENTERO-Screen

Kat. č.: MLT00009

Pro mikrobiologii

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Cat. No.: MLT00009

For microbiology

The kit ENTERO-Screen is intended for the identification of enterobacteria by means of eight biochemical tests: glucose (anaerobic), acetoin, phenylalanine, indol, sucrose, urease, lysine, or-nithin. The composition of tests and the 4 hours incubation period determine the kit for usage as Salmonella screening test in food industry and for identification of enterobacteria mainly in urinary samples. However it can be used for other clinical material as well. The kit enables to perform 36 determinations. The identification can be supplemented by tests for cytochrome oxidase (OXIttest), *E. coli* (COLIttest) and pyrrolidonylamidase (PYRAtest).

The kit ENTERO-Screen contains:

- 3 microtitration plates (for identification of 12 strains each) with desiccant
- Instructions for use including the differentiation table and index
- Colour scale for ENTERO-Screen kit
- 3 polyethylene bags for incubation
- Storage bag (for storage of an open plate), 1 pc
- 36 record sheets
- Lid

Storage, expiration:

ENTERO-Screen should be stored at (+2 to +25)°C. Expiration is indicated in each package.

Recommended working procedure for ENTERO-Screen

Required material for performing the identification (not included in the kit):

- Reagent for ACETOIN test (Cat. No. MLT00016 – 270 determinations)
- Reagent for PHENYLALANINE test (Cat. No. MLT00017 – 230 determinations)
- Reagent for INDOLE test (Cat. No. MLT00020 – 310 determinations)
- Paraffin oil sterilized (Cat. No. MLT00042 – 120 determinations)
- Petri dishes with the cultivation medium
- Test tubes (100x15) mm with 1 ml of a sterile physiological saline
- Instrument DENSILAMETER II (Cat. No. INS00062)
- Automatic micropipette 0.1 ml, sterile tips.
- Thermostat (35–37)°C.
- Usual equipment of a microbiological laboratory (loops, markers, burner).

For supplementary test performing (not included in the kit):

- OXIttest (Cat. No. MLT00039 – 50 determinations)
- COLIttest (Cat. No. MLT00035 – 50 determinations)
- Reagent for INDOLE test (Cat. No. MLT00020 – 120 determinations)
- PYRAtest (Cat. No. MLT00040 – 50 determinations)
- Reagent for PYRAtest (Cat. No. MLT00023)

For results evaluation (not included in the kit):

- Code Book for ENTERO-Screen - located at www.erbalachema.com
- The ErbaExpert Identification Program

Caution:

- For professional use only

Follow principles for working with infectious material!

Isolation of cultures:

- Perform the isolation of cultures by usual technique on a suitable cultivation medium (MacConkey agar, blood agar etc.).
- Perform Gram staining and the oxidase test (by using the OXIttest).

Note:

- In the case of using the ENTERO-Screen for the uropathogens identification, it is recommended to complement the procedure by performance the COLIttest, resp. PYRAtest.
- Procedure for the COLIttest and PYRAtest is described in the working instructions of the individual tests.

Preparation of inoculum:

Make a bacterial suspension equal to McFarland No. 3 turbidity scale in physiological saline from pure 24 hours old culture.

Culture purity control:

If required to confirm purity of the suspension or growth properties of the culture, streak-out a sample from the inoculated physiological saline on non-selective medium and incubate at 37 °C for 24 hours.

Preparation of the ENTERO-Screen plate:

- Open the aluminium sachet close to the weld and take out the plate.
- Cut off required number of triple strips from plate.
- Remove the adhesive tape from individual strips and insert them into the prepared frame. In case you work with MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely.

- Record number of the strains or isolates to be examined on the appropriate strips
- Put the rest of the plate with desiccant to the storage bag and store it in a refrigerator for further use; protect it from humidity. It is recommended to spend the rest of plate till 4 weeks after the first use.

Note:

Any uneven distribution of substrate in the well does not affect the functionality of the test.

Inoculation:

- Shake the bacterial suspension thoroughly before use.
- Inoculate 0.1 ml of the suspension into all microwells of the respective strip.
- Avoid the contamination of adjacent microwells during the inoculation process.
- After the inoculation, drop the paraffin oil into:
 - the wells in the column H (test glucose) – 4 drops of the paraffin oil
 - the wells in the columns E, C, B, A (tests indole, urease, lysine, ornithine) – 2 drops of the paraffin oil.

Note:

- Inoculate control strains simultaneously with any series of unknown strains as well as when using a new batch of the plates ENTERO-Screen to control the colour expression of positive and negative reactions.

Note:

Clean the inside of the lid by ethanol just before the use.

Incubation:

- Insert the frame with inoculated strips into the polyethylene bag.
- Fold the open end of the bag under the plate to prevent evaporation.
- Incubate at (+35 to +37 °C) for 4 hours.

Reading reactions:

- Before reading drop reagents to the following tests:
 - acetoin (VPT, wells in the column G) – one drop of the reagent for VPT I and VPT II respectively
 - indole (IND, wells in the column E) – two drops of the reagent for IND.
- Incubate the plate for other 30 minutes at (35–37) °C
Then drop the reagent to the phenylalanine test (PHE, wells in the column F) – one drop of the reagent for PHE and read this reaction immediately (the colouration of positive reaction disappears within 2–3 minutes).
- Read all remaining tests and record the results in the record sheet.

Note:

- Use the table "Interpretation of reactions", the Colour scale for ENTERO-Screen or proceed in accordance with colour reactions of the control strains to evaluate the colour reactions properly.

Results interpretation:

Screening of salmonellas

- ENTERO-Screen plate is designed to eliminate quickly bacteria that do not belong to the *Salmonella* genus:
 - the negative test GLU, resp. positive reaction with some of the tests VPT, PHE, IND, SUC and URE indicate that the isolate does not belong to salmonellas.

Note:

- Certain cultures of salmonellas may give a negative or delayed reaction LYS or ORN, in the case of suspicious identification use ENTEROtest.

Identification of enterobacteria:

- To evaluate the results use Differentiation table or Index included in the working instructions.
- To complete identification evaluate the culture as a complex, considering the kind of colonies, the microscopy, the lactose reaction, origin of the isolate etc.
- The culture that could not be identified successfully by ENTERO-Screen should be identified by the kits ENTEROtest.

Note:

- In the case of a culture with the negative test "Glucose in the anaerobic method", it is highly probable that it is a non-fermenting bacterium. For identification of the gram-negative non-fermenting bacteria, the kit NEFERMtest 24 is recommended.

Disposal of the used material:

- Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.
- Put paper packaging waste to recycling.

The most frequent causes of identification failure:

- A mixed or contaminated culture.
- Inoculum of a low density or a small volume.
- Inoculum has contaminated adjacent strips.
- The respective tests have not been overlayed by the paraffin oil.
- The reagent was dropped into the adjacent well.
- Failure to follow the recommended working procedure.

In case of an atypical strain or representative of a species that has not been included in the Differentiation table and Index, perform the identification by ENTEROtest kits.

Identification by means of Index:

- The Index includes typical results of the mostly often isolated enterobacteria from urinary infections.
- Both the results and interpretations are arranged in the Index according to the profile values.

Profile calculation:

- The results in record sheet are subdivided into three groups.
- In each group, the tests are marked with numbers 1, 2 and 4.
- When calculating the profile, sum up the values of numbers with tests delivering positive reactions in any group of tests (whereas the tests giving negative results are marked with zero).
- Record the sum into the respective field of the row "Profile" in the record sheet.
- The resulting number consisting of three digits = Profile.

Example of the Profile calculation:

H GLU 1	G VPT 2	F PHE 4	E IND 1	D SUC 2	C URE 4	B LYS 1	A ORN 2	
+	+	-	-	+	+	+	-	- ENTERO-Screen
1	2	(0)	(0)	2	4	1	(0)	- numerical test values
								- reaction results
								- sum of values
								- profile

Performance:

The kit was tested on a set of 120 strains.
The identification of 93 % strains was correct.
7 % of the strains were not identified.

Quality control of tests:

The quality control of the kits is performed systematically at various stages of their production. The batches are checked by tests on standard bacterial cultures. For those who wish to perform their own quality control tests, cultures mentioned in the table **Control strains** are recommended. **These strains are used to check the functionality of the kit, not to check the accuracy or success of the identification!**

Enterobacter cloacae CCM 1903 (ATCC 10699)

Proteus proteus sp. CCM 1799

Edwardsiella tarda CCM 2238 (ATCC 15947)

Control strains

CCM No.	Column/Abbreviation of the test							
	H GLU	G VPT	F PHE	E IND	D SUC	C URE	B LYS	A ORN
1903	+	+	-	-	+	-	-	+
1799	+	-	+	+	+	+	-	-
2238	+	-	-	+	-	-	+	+

Legend: + = positive reaction - = negative reaction

The strains are supplied in freeze-dried ampoules by the CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz.

Caution:

It is necessary to use fresh isolates of the CCM strains each time when kit checked up for functionality. These strains are intended for check-up of the functionality of the kit, not for check-up of accuracy or effect of the identification.

Health protection:

Components of the kit are not classified as dangerous.

USED SYMBOLS

Catalogue number

In vitro diagnostics

Manufacturer

See instruction for use

Lot number

Storage temperature

Expiry date

Interpretation of reactions

Column	Test	Code	Reaction	
			positive	negative
H	Glucose	GLU	yellow, yellow-to-green	green
G	Acetoin	VPT	red, pink	colourless
F	Phenylalanine	PHE	dark green	yellow
E	Indole	IND	red, pink	yellow
D	Saccharose	SUC	yellow, yellow-to-green	green
C	Urease	URE	red, red-to-orange	yellow, yellow-to-orange
B	Lysine	LYS	blue, blue-to-green	yellow-to-green, green
A	Ornithine	ORN	blue, blue-to-green	yellow-to-green, green

Diferentiation table

v.1.0

ENTERO-Screen								Identification	Additional tests			
H GLU	G VPT	F PHE	E IND	D SUC	C URE	B LYS	A ORN		GLR	PYR	LAC	ONP
+	+	-	+	+	d	+	-	Klebsiella oxytoca	-	+	+	+
+	+	-	-	+	-	+	+	Enterobacter aerogenes	-	+	+	+
+	+	-	-	+	(-)	+	+	Serratia marcescens	-	+	-	+
+	+	-	-	+	(+)	+	-	Klebsiella pneumoniae	-	+	+	+
+	+	-	-	+	d	-	+	Enterobacter cloacae	-	d	+	+
+	d	+	-	(-)	+	-	+	Proteus mirabilis	-	-	-	-
+	-	+	+	-	+	-	+	Morganella morganii	-	-	-	-
+	-	+	+	+	+	-	-	Proteus vulgaris	-	-	-	-
+	-	+	+	(-)	+	-	-	Providencia rettgeri	-	-	-	-
+	-	+	+	d	d	-	-	Providencia stuartii	-	-	-	-
+	-	-	+	d	(+)	-	+	Citrobacter koseri	-	+	d	+
+	-	-	+	(-)	(+)	-	+	Citrobacter amalonaticus	-	+	d	+
+	-	-	+	d	-	(+)	d	Escherichia coli	(+)	-	+	+
+	-	-	-	d	d	-	(-)	Citrobacter freundii komplex	-	+	d	+
+	-	-	-	-	-	+	+	Salmonella sp. ¹⁾				

Legend: + = positive reaction - = negative reaction d = variable reaction

(+) = mostly positive reaction (-) = mostly negative reaction

¹⁾: Certain cultures of salmonellas may give a negative or a delayed reaction LYS or ORN

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