NEISSERIAtest

Kat. č.: MLT00011

Pro mikrobiologii

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Cat. No.: MLT00011

For microbiology

The NEISSERIAtest is intended for the identification of bacteria of the genus *Neisseria*, especially *N. gonorrhoeae*, *N. meningitidis* and *Moraxella (Branhamella) catarrhalis*. The kit allows you to perform identification of 36 bacterial strains, using 7 biochemical tests. Identification is complemented with tests for detection of oxidase and β -galactosidase, supplied in the form of dg. strips OXItest and ONPtest. The use of divided plates makes it possible to use only the necessary part of the plate, corresponding to the number of tested strains.

- The kit NEISSERIAtest contains:**
- 3 microtitration plates (for the identification of 12 strains each) with desiccant
 - 3 polyethylene bags for incubation
 - Storage bag (for storage of an open plate), 1 pc
 - Instructions for use including the differentiation table,
 - Colour scale for NEISSERIAtest
 - 36 record sheets
 - Lid

Storage, expiration:

The NEISSERIAtest kit should be stored at (+2 to +8) °C. Store away from direct sunlight. The expiration date is indicated on the package.

Recommended procedure

Material required to perform a test (not included in the NEISSERIAtest kit):

- Suspension medium for NEISSERIAtest (Cat. No. MLT00025 – 18 determinations)
- Lugol iodine solution
- Petri dishes with cultivation medium suitable for neisseria
- Incubator, equipment for cultivation in a CO₂ enriched atmosphere
- Automatic micropipette (set to 0.1 ml) plus sterile tips
- Regular microbiological laboratory equipment (loops, marker, burner, scalpel)
- McFarland densitometer

Material required to perform a supplementary test (not included in the NEISSERIAtest kit):

- OXItest strips (Cat. No. MLT00039 – 50 determinations)
- Reagent for OXIDASE test (Cat. No. MLT00022 – 250 determinations)
- ONPtest strips (Cat. No. MLT00038 – 50 determinations)
- V+K DISK (Cat. No. MLT00087 – 100 determinations)

For results evaluation (not included in the kit):

- The ErbaExpert Identification Program

Caution:

For professional use only

Follow the principles for working with infectious material!

Isolation of cultures:

Isolation of cultures is performed on chocolate Mueller-Hinton agar, blood agar No. 3 with growth enrichment, or on agars with an inhibitory supplement (Thayer-Martin agar), or on other agars recommended for *Neisseria* cultivation.

- Heat the agar dishes in the thermostat to a temperature of (35-37) °C before inoculation
- Incubate the plates with the inoculated culture in an environment with an increased CO₂ content and increased humidity, at a temperature of (35-37) °C for 18-24 h
- From fresh grown culture perform a Gram stain, oxidase test, or ev. catalase test, microscopy
- Gram-negative cocci ev. diplococci with positive oxidase and catalase identify using the kit NEISSERIA test (*N. elongata* forms rods, catalase negative)

In the case of working with material in which mixed cultures may occur (e.g. from respiration infection, when using a non-selective medium), we recommend:

- inoculate a well-isolated colony from the primary culture
- perform its repeated cultivation before self-identification using NEISSERIAtest

For the selective isolation of *Neisseria meningitidis* from the throat diagnostic disc V+K DISK, Cat. No. MLT00087, can be used. In the zone of antibiotic diffusion grow only *N. meningitidis*, *N. lactamica* and yeasts. Their differentiation:

- yeast - by microscopy
- *N. lactamica* - using the ONPtest, Cat. No. MLT00038 - up to 1 h ONPG gives a positive reaction (yellow color), ONPG negative strains can be identified using the NEISSERIAtest
- *N. meningitidis* (ONPG negative) - identification using the NEISSERIAtest with results within 4 hours

Preparation of the NEISSERIAtest plate:

- Open the aluminum sachet close to the weld and remove the plate
 - Using a scalpel, cut the appropriate number of rows (strips) of the plate, corresponding to the number of tested strains (1 strip, i.e. 1x8 tests, to identify one strain). In case you are working with the MIKROLATEST® kit for the first time and an empty frame is not available, use the frame of the first plate. The unused strips of the first plate put into the storage bag freely
 - Remove the cut strips from the panel, remove the protective Al foil and place the strips in the prepared empty frame
 - Record the numbers of the tested strains or isolates on the appropriate strips
 - Put the rest of the unused plate with the desiccant to the bag for storage of open plate, enclosed with the kit, and store in the refrigerator for further use; make sure that the plate is protected against humidity. It is recommended to spend the rest of plate till 4 weeks after the first use
- Note:** Any uneven distribution of the substrate in the well does not affect the functionality of the test

Preparation of inoculum:

- Preheat the suspension medium and the NEISSERIAtest kit to (35–37) °C
- From the culture grown on the agar medium, prepare a suspension in the suspension medium with a density corresponding to at least level 3 of the McFarland turbidity scale using an inoculation loop or possibly a cotton swab, thoroughly homogenize the suspension
- **Use a 24 h culture of the tested strain;** the use of an older culture may lead to false negative results
- At the same time, make a cross smear from the suspension on the agar plate to check the purity, possibly for further culture testing
- In the case of using physiological solution for the ONPtest, prepare a suspension in the physiological solution of the same density as for the inoculation of the NEISSERIAtest

Inoculation:

- Inoculate the strips as soon as possible; not later than 15–30 min after preparation of the suspension
 - Homogenize the suspension thoroughly before use
 - Inoculate 0.1 ml of the suspension into each well of the respective strip
 - Do not dispense the inoculum too vigorously to prevent aerosol contamination of the wells of adjacent strip
 - Cover the inoculated strips with lid
 - To the rest of the suspension in the suspension medium, possibly into the suspension in physiological solution, insert a strip for the detection of β -galactosidase (ONPtest) using annealed tweezers
- Note:** If you use the lid to cover the plate during work, wipe its inside with ethanol before use.
- With each series of unknown strains and always when using a new batch of NEISSERIAtest plates, simultaneously inoculate control strains to verify the color expression of the tests

Incubation:

- Insert the NEISSERIAtest plate with inoculated strains into the incubation PE bag
- Fold the open end of the bag under the plate to prevent the inoculum from drying out
- Leave the tubes with the ONPtest in the rack
- Incubate the NEISSERIAtest and the ONPtest in the rest of the suspension at (35–37) °C for 24 hours in aerobic conditions, in the case of using a sufficiently thick suspension of the strain, it is possible to perform preliminary assessment after 4 hours of incubation (without SPS test, well A; read SPS only after 24 hours of incubation)

Reading:

- Record the results of reactions in wells G - B and ONPtest into record sheet
- Compare the results of sugar utilization assay (wells G - D) with negative control (well H)
- **Any change of the colour** reactions of these sugars in comparison to negative control in well H means positive reaction
- To evaluate the colour reactions of the tests follow the table "Interpretation of reactions" and/or the colour reaction of the control strains
- Add one drop of the Lugol reagent into well A (test SPS) and record the result of the reaction

Identification:

- For the identification refer to the "Differentiation table".
- To complete the identification take into consideration all the results available, i. e. the source of culture, appearance and consistency of colonies, growth on selective media etc.
- If you have failed to identify the culture repeat the procedure as above.

Disposal of used material:

- After use, place the inoculated strips of the plate in an infectious material container and autoclave or destroy by incineration
- Put paper packaging waste to recycling

The most frequent causes of identification failure:

- Mixed or contaminated culture
- Using inoculum of low density or small volume of inoculum
- Inoculum has contaminated adjacent strip
- Failure to follow the recommended procedure
- There may be a species or strain whose data are not included in the "Differentiation table"

Performance:

- The kit was tested on a set of 112 strains.
- The identification of 98% strains was correct.
- 2% of the strains were not differentiated.

Quality control of NEISSERIAtest:

The quality of the chemicals used for the production of NEISSERIAtest plates is verified by a standard test procedure. The produced series of plates are also checked by a functional test using control bacterial strains. To work with NEISSERIAtest plates at your workplace, we recommend using the control strains listed in the table (see below). Also for routine diagnostics, we recommend using these standard test strains to verify the correctness of the methodical procedure, the course of the tests and the color expression of the reactions. **Control strains** can be recommended to be used with each series of unknown strains and always when using a new batch of the kit, respectively according to the laboratory's validation procedure. To check the functionality of the kit, it is always necessary to use fresh isolates of the control strains. **Attention - these strains are only used to check the functionality of the kit, not to check the correctness or success of the identification!**

Control strains

| Strain | CCM | Well. No. | | | | | | | | Strip ONPtest |
|-----------------------------------|------|-----------|----------|----------|----------|----------|----------|----------|----------|---------------|
| | | H NEC | G GLU | F MLT | E FRU | D SUC | C GGT | B TRB | A SPS | |
| <i>Neisseria subflava</i> | 3482 | – | + | + | + | + | + | – | + | – |
| <i>Neisseria lactamica</i> | 4392 | – | + | + | – | – | – | – | – | + |
| <i>Moraxella (B.) catarrhalis</i> | 4391 | – | – | – | – | – | – | + | – | – |

Explanations: +..... positive reactions
 –..... negative reactions
 NEC..... negative control

These strains are supplied in freeze-dried ampoules by the CCM – Czech Collection of Micro-organisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, tel. 549 491 430, fax 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz

References:

Morse, S. A., and Knapp, J.S.: The genus *Neisseria*. In: The Prokaryotes. 2nd. ed. Edited by E. Ballows, H. G. Trüpper, M. Dworkin, W. Harder, and K. H. Schleifer. Vol. 3. p. 2495–2529. Springer Verlag, New York, Berlin, Heidelberg 1991.

Health protection:

Components of the kit do not contain any dangerous substances.

Interpretation of reaction

| Well No. | Test | Code | Reaction | |
|----------|--------------------------|------|--------------------------------|-------------------------------|
| | | | Positive | Negative |
| H | Negative control | NEC | – | red to orange-red |
| G | Glucose | GLU | yellow or pale orange | red to orange-red |
| F | Maltose | MLT | yellow or pale orange | red to orange-red |
| E | Fructose | FRU | yellow or pale orange | red to orange-red |
| D | Sucrose | SUC | yellow or pale orange | red to orange-red |
| C | λ-glutamyl transferase | GGT | yellow or pale yellow | colourless |
| B | Tributyrin | TRB | yellow or yellow-orange | red or orange-red |
| A | Polysaccharide synthesis | SPS | black, dark blue or dark brown | yellow, brown or brown-yellow |
| Strip | | | | |
| OXItest | cytochrome oxidase | OXI | dark blue or blue | grey |
| ONPtest | β-galactosidase | ONP | yellow or pale yellow | colourless |

Differentiation table

| NEISSERIAtest | | | | | | | | | Supplementary tests | | |
|----------------------------|-----|-----|-----|-----|-----|-----|-----|-----|---------------------|-----|--------------------------------|
| Well No.: | H | G | F | E | D | C | B | A | | | Growth on Thayer-Martin medium |
| Test: | NEC | GLU | MLT | FRU | SUC | GGT | TRB | SPS | ONP | DNA | |
| <i>N. gonorrhoeae</i> | – | + | – | – | – | – | – | – | – | – | + |
| <i>N. meningitidis</i> | – | + | + | – | – | + | – | – | – | – | + |
| <i>N. lactamica</i> | – | + | + | – | – | – | – | – | + | – | + |
| <i>N. polysaccharea</i> | – | + | + | – | d | – | – | + | – | – | + |
| <i>N. sicca</i> + | – | + | + | + | + | d | – | + | – | – | – |
| <i>N. mucosa</i> + | | | | | | | | | | | |
| <i>N. subflava</i> | – | + | + | d | d | d | – | d | – | – | – |
| <i>N. flavescens</i> | – | – | – | – | – | – | – | + | – | – | – |
| <i>N. cinerea</i> ++ | – | – | – | – | – | – | – | – | – | – | – |
| <i>N. elongata</i> ++ | | | | | | | | | | | |
| <i>M. (B.) catarrhalis</i> | – | – | – | – | – | – | + | – | – | + | – |

Explanations: + = 75–100 % of positive reactions
 – = 0–25 % of positive reactions
 d = 26–74 % of positive reactions
 NEC = negative control
 DNA = deoxyribonuclease

Notes: + *N. sicca* and *N. mucosa* are not differentiated by NEISSERIAtest
 ++ *N. cinerea* cells occur as plump cocci arranged in pairs or in scattered clusters
N. elongata cells occur as short and slender rods often arranged as diplobacilli or in short chains

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