mOXItest

Kat. č.: MLT00030

Pro mikrobiologii

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Cat. No.: MLT00030

For microbiology

The mOXItest is intended for the detection of the cytochrome oxidase activity, for the purpose of the microbiological control of samples applying the membrane filter method. The preparation can be used for the control of drinking water, within the food industry etc. One package of the mOXItest enabling 15 determinations.

Principle:

The membrane filter with bacterial colonies shall be transported on the mOXItest disc surface, the activity of the cytochrome oxidase being detected by the colour reaction of N,N-dimethyl-1,4-phenylenediamine with α -naphthol which form indophenol blue.

A package of mOXItest contains:

- 5 bags with discs for the oxidase detection and with a siccative, each bag containing 3 pcs of discs for 3 determinations
- Working leaflet

Caution:

For professional use only.

Storage:

mOXItest should be stored at the temperature of (+2 to +8) °C. Store away from direct sunlight. The expiry date is indicated on each package.

Required material:

- Pipette 0.6 ml, distilled water
- Empty Petri dishes with a diameter of 6 cm
- Tweezers
- Usual equipment of a microbiological laboratory (burner, markers etc.)

Recommended procedure:

- Insert 1 disc of mOXItest into an empty Petri dish.
- Moisten the disc with 0.6 ml of distilled water (sterile provided, you would like to restreak the bacteria colonies tested for the next cultivation).
- Applying the tweezers, carry over carefully the membrane filter with the colonies on the disc so that there are no air bubbles left between the supporting plate and the filter there (it is recommended to lay the filter upon the disc gradually, from one filter edge to the other one).
- After approx. 1–2 minutes evaluate visually the oxidase reaction with individual colonies according to the data delivered in the table: **Interpretation of the reaction**. In the case of evaluating the reaction after a period exceeding 2 minutes, the false positive reaction may occur.

Note:

- When working with mOXItest, avoid any contact of your hands with the disc surface, after having finished the work, wash up your hands thoroughly using water and soap.
- Insert the used material into the vessel intended for the infectious material and autoclave or destroy it by incineration.
- In the case of not using up the bag contents immediately after opening, seal-in the discs left together with the siccative back into the bag, resp. store them dry and kept in the dark, nevertheless not longer than 1 week after having the bag open.

Interpretation of the reaction:

Reaction	Coloration delivered
positive	blue, pale blue coloration
negative	no blue coloration

Standard control strains:

For to check and control the quality of the mOXItest, the control strains given in the following table can be applied: using the inoculation loops, apply the growth of the control strain from the agar medium directly to the mOXItest disc surface. Read the colour reaction after 1–2 minutes of incubation.

Control strains:

Strain	Reaction
<i>Pseudomonas aeruginosa</i> CCM 1960	positive
<i>Enterobacter sakazakii</i> CCM 3461	negative

These strains are supplied in freeze-dried ampoules the Czech Collection of Microorganisms,

Masaryk University, Kamenice 5, 602 00 Brno, Czech Republic. ccm@sci.muni.cz.

Phone: (+420) 549 491 430, Fax: (+420) 543 247 339, e-mail: ccm@sci.muni.cz.

The strains are delivered in a lyophilized form or on gelatinized discs.

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