COLItest

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Dátum revízie: 13. 8. 2020

COLItest

Cat. No.: MLT00035
The identification strips COLItest are intended for quick identification of Escherichia coli bacteria, isolated from clinical materials, foodstuffs, waters etc. The identification is based on the determination of β-glucuronidase activity and indole formation. The reaction can be after only 4 hours of incubation.

Principle:
The strip of COLItest is inserted into the suspension of the tested strain in saline solution and incubated at 37 °C. Reaction substances are released into the suspension where the reactions are read after time of incubation. The β - glucuronidase enzyme splits 4 - methylumbelliferyl - β - D - glucuronide; released 4 - methylumbelliferon gives blue fluorescence under an UV lamp. Production of indole from L-tryptophan is detected by red colouring after addition of indole specific reagent (solution of p-dimethylaminobenzaldehyde).

A package of the COLItest contains:
• 50 test strips for 50 determinations
• working leaflet

Storage:
In refrigerator at (+2 to +8) °C. The time of expiration is marked on each package.

Caution:
• For professional use only

Materials:
• COLItest strip
• Sterile saline solution 0.5–1 ml in test tubes
• Incubator 37 °C
• Long wavelength UV lamp, about 360 nm (for example commercially available battery – operated UV lamps)
• Reagent for INDOL test cat. No. MLT00020
• Instrument DENSILAMETER II cat. No. INS00062

Recommended procedure:
• Use COLItest for coliform bacteria and/or gramnegative, oxidase negative rods
• Prepare suspension at 0.5–1 ml of sterile saline solution from the tested strain
• The density if the suspension should be equal to number 3 McFarland's barium sulphate standard
• Insert one COLItest strip into the tube with suspension - the paper zone of the strip has to be soaked into the suspension
• Incubate at 37 °C
• After 4 hours of incubation:
- first, read β-glucuronidase reaction (GLR) and record the result
- use UV lamp 360 nm, read the GLR test in a dark place; the positive reaction gives blue fluorescence
- after GLR test, read indole test (IND); add 4–5 drops of Reagent for IND the positive reaction is indicated by red colouring

Note:
• The GLR test can be preliminary read after 1 hour of incubation

Waste disposal:
Used strip should be treated as potentially infectious and should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging waste to recycling or in municipal waste.

Interpretation of the reactions			Explanations: GLR - β-glucuronidase test IND - indole production test
Reaction	GLR (fluorescence)	IND (colour reaction)	
positive	blue fluorescence	red to pink colour	
negative	no fluorescence	yellow to yellowish colour	

Recommendations:
• The glass tubes are free of any fluorescence test the tubes under UV lamp before starting the test
• Always, use negative control; insert 1 strip of COLItest into the pure saline solution and incubate it together with the tested strains suspensions; the negative control does not exhibit any fluorescence
• It is possible to use recommended control strains - see table "Control strains"
Evaluation:
Typical Escherichia coli strains gives GLR + IND + reactions. Some strains of Escherichia coli can give GLR+ IND - or GLR- IND+ reactions. Confirm these strains by further identification. Combination of GLR-IND+ is typical for E. coli O157 H:7.
Note:
• About 95 % strains of Escherichia coli are β-glucuronidase positive; positive indole reaction increases sensitivity and specificity of E. coli identification
• The COLItest is also suitable for differentiation of other groups of bacteria, for example for detection of β-glucuronidase activity within staphylococci etc.

Strain	CCM	ATCC	GLR	IND
Escherichia coli	3988	10536	+	+
Serratia marcescens	303	13880	-	-

References:
Chang, G. W. et al. Appl. Environ. Microbiol. 55, 335-339, 1989
Dealler, S. F. et al. Med. Lab. Sci. 49, 12-15, 1992
Dominguez, A. et al. Diagn. Microbiol. Infect. Dis. 15, 291-294, 1993
Edberg, S. C., and Kotnick, C. M. J. Clin. Microbiol. 24, 368-371, 1986
Kämpfer, P. et al J. Clin. Microbiol. 29, 2877-2879, 1991,
Manafi, M. Zbl. Hyg. 189, 225-234, 1989
Manafi, M. and Kneifel, W. Acta Microbiol. Hung. 38, 293-304, 1991
Rice, E.W. et al. Appl. Environ. Microbiol. 56, 1203-1205, 1990

The strains are delivered in a lyophilized form or on gelatinized disks.

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