ONPtest

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Cat. No.: MLT00038

For microbiology

Principle:
The ONPtest is intended for the detection of bacterial β-galactosidase. A package of the ONPtest which contains 50 test strips enables 50 determinations. The β-galactosidase enzyme splits o-nitrophenyl-β-D-galactoside; released o-nitrophenol gives yellow colour, indicating thus positive reaction.

A package of the ONPtest contains:
• 50 test strips for 50 determinations
• working leaflet

Storage:
In refrigerator at (+2 to +8) °C. The time of expiration is indicated on each package.

Caution:
• For professional use only

Recommended procedure:
• Use the pure culture of the tested strain
• Take a tube (100 x 15 mm) with 0.6 ml of physiological saline
• Emulsify the pure culture with a loop in the tube by rubbing against the side of the tube and mixing
• The density of the suspension should be equal to number 2 of McFarland's turbidity scale
• Insert the strip into the tube – the zone of the strip has to be soaked into the suspension
• Incubate
– at 37 °C for 24 hours for Enterobacteriaceae or Neisseria
– at 37 °C for 24 - 48 hours for Gram-negative non-fermenting bacteria
• Read the reaction visually after time of incubation – evaluate the colour of the suspension

Waste disposal:
Used strip should be treated as potentially infectious and should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging waste to recycling or in municipal waste.

Note:
The ONPtest can be used as additional test for MIKROLATEST® kits; follow instructions of the working leaflet in the kit.

Interpretation of the reaction

Test	Code	Reaction	
β-galactosidase (ONPG)	ONP	positive	negative
		yellow, light yellow	colourless

It is recommended to use standard control strains for checking of the reaction and interpretation of the colours of the reaction.

Standard control strains:

Strain	CCM	ATCC	ONP
Enterobacter sakazaki	3461		+
Pseudomonas aeruginosa	1960	10145	-

These strains are distributed by
Czech Collection of Microorganisms, Masaryk University,
Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ
Tel: +420 549 491 430, Fax: +420 549 498 289 http://www.sci.muni.cz/ccm, e-mail: ccm@sci.muni.cz

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