MIC NEFERM

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M-MIC/PI/04/25/K/INT

**Cat. N.: MLT00046****For microbiology**

The kit is designed to test antimicrobial susceptibility of Gram-negative non-fermentative bacteria on the basis of MIC (minimal inhibitory concentration) determination, i.e. the lowest concentration, which inhibits bacterial growth. The kit contains 10 examinations (plates). The test is based on the rehydration of antibiotics in the wells with Muller Hinton II broth and addition of bacterial suspension. The results are read visually or by reader after 16–20 hours of incubation.

The kit contains:

- 10 plates for examination
- A lid (non-sterile)
- 10 pc of PE bags

Storage and expiration of the kit:

It is recommended to store the kit at (+2 to +25) °C. The date of expiration is indicated on each package. Leave plate at room temperature at least for 30 minutes before you open it to avoid water condensation. After the aluminium package is opened, don't leave opened plates unprotected!!! Exposition to air humidity leads to antibiotic activity failure!!!

Material required to perform a test, not included in the kit:

- Sterile physiological solution (unbuffered)
- Muller Hinton broth II cations adjusted (e.g. suspension medium MIC Cat. N. MLT00070)
- Ethanol
- Sterile tubes
- Inoculator ErbaDip Cat. N. 50004456
- Sterile Petri dishes
- Sterile basins 60 ml Cat. N. 50004457
- A stepper or multichannel pipette for dosage of 100 µl
- A pipette for dosage of 60–100 µl
- Densitometer (e.g. DENSILAMETER II Cat. N. INS00062)
- Incubator 35±2 °C
- Regular microbiological laboratory equipment (loops, marker, burner, etc.)
- Program ErbaExpert

Caution: The kit is for professional use only! Respect the rules for work with infectious material!

Instructions for Use**Preparation of bacterial suspension and inoculation (recommended procedure):****A) Inoculation with inoculator**

- 1) Remove a plate from aluminium bag and remove aluminium cover. Mark the frame with a type of kit (NEFERM) to avoid mistake in reading results after incubation. Record number of examined strain on the plate. Fill 100 µl of suspension medium MIC into each well.
- 2) Prepare a tube with 12 ml of physiological solution. Add 100 µl of suspension medium MIC to decrease surface tension.
- 3) Remove few colonies from 18–24 hour culture on blood agar and prepare a bacterial suspension of density of 0.5 on McFarland scale in physiological solution. It is necessary to give extra attention to achieving a smooth suspension of 0.5 McFarland when preparing inoculum of *Pseudomonas* spp. The growth method can be used alternatively: 3–5 colonies are transferred into a tube with suitable broth medium and incubated until the turbidity of 0.5 McFarland is achieved or exceeded. The turbidity of the broth culture is adjusted with sterile saline or broth to 0.5 McFarland standard.
- 4) Pour the bacterial suspension into a sterile Petri dish.
- 5) Use sterile inoculator to inoculate the plate: dip inoculator into Petri dish with ethanol and flame it. Dip the cooled inoculator into a Petri dish with prepared bacterial suspension. A thin film of bacterial suspension is adhered to metal spikes of inoculator. Transfer inoculum to the first half of the plate by dipping into wells and careful mixing. Make a new dip into the Petri dish with prepared bacterial inoculum and inoculate the second half of the plate.

B) Inoculation with pipette

- 1) Prepare a tube with 2 ml of physiological solution.
- 2) Remove few colonies from 18–24 hour culture on blood agar and prepare a bacterial suspension of density of 0.5 on McFarland scale in physiological solution. It is necessary to give extra attention to achieving a smooth suspension of 0.5 McFarland when preparing inoculum of *Pseudomonas* spp. The growth method can be used alternatively: 3–5 colonies are transferred into a tube with suitable broth medium and incubated until the turbidity of 0.5 McFarland is achieved or exceeded. The turbidity of the broth culture is adjusted with sterile saline or broth to 0.5 McFarland standard.
- 3) Place 60 µl of bacterial suspension into a tube with 13 ml of suspension medium MIC, homogenise well.
- 4) Remove a plate from aluminium bag and remove aluminium cover from the plate. Mark the frame with a type of kit (NEFERM) to avoid mistake in reading results after incubation. Record number of the examined strain on the plate.
- 5) Inoculate each well of the plate with 100 µl of bacterial suspension prepared in suspension medium MIC.

Note: Process a plate within 60 minutes after removing it from aluminium bag.

Incubation:

Insert the inoculated plate into a PE bag. Fold the open end of the bag under the plate to prevent evaporation during the incubation. Incubate the plate at 35±2 °C for 16–20 hours. If MIC can not be evaluated in slowly growing strains (e.g. *Stenotrophomonas maltophilia*, *Burkholderia cepacia* and *Pseudomonas aeruginosa* from cystic fibrosis) extend the incubation for another 20 hours at room temperature.

Evaluation:

Remove the plate from the PE bag. To read the growth in the microwells, choose a way which is the most convenient for you:

- 1) Read against a grey background or against plate layout in instructions.
- 2) Read against natural or artificial dispersed light.
- 3) Usage of magnifying glass is not recommended.
- 4) Evaluate test using reader and identification software.

Please read with attention!

You must see a growth in the control well (K)! If the growth is not present, the test MUST NOT be evaluated! The MIC is the lowest concentration of antibiotic in a well where no visible growth of the organism is observed. Exception: With Trimehtoprim/sulfamethoxazol, a well with ≥ 80 % growth inhibition compared to the growth control is considered as MIC. In Colistin may occur very weak growth, which is also considered a positive result. In case of doubt about the presence of growth in the wells, use ErbaScan for reading. Beware to differentiate grains of growth from media bubbles. Record the results.

Tab. 1: Plate layout: antibiotics dilution series (in mg/l)

	1 AMS	2 PIP	3 PIT	4 CAZ	5 AZT	6 MER	7 GEN	8 AMK	9 COL	10 CIP	11 TGC	12 T/S
A	128/64	128	128/4	16	16	16	32	64	16	8	8	4/76
B	64/32	64	64/4	8	8	8	16	32	8	4	4	2/38
C	32/16	32	32/4	4	4	4	8	16	4	2	2	1/19
D	16/8	16	16/4	2	2	2	4	8	2	1	1	0.5/9.5
E	8/4	8	8/4	1	1	1	2	4	1	0.5	0.5	0.25/4.75
F	4/2	4	4/4	0.5	0.5	0.5	1	2	0.5	0.25	0.25	0.12/2.38
G	2/1	2	2/4	0.25	0.25	0.25	0.5	1	0.25	0.12	0.12	0.06/1.19
H	1/0.5	1	1/4	0.12	0.12	0.12	0.25	0.5	0.12	0.06	0.06	K

Interpretation:
Depending on national or laboratory standards, current interpretation tables by EUCAST (1), CLSI M100 (2) or other interpretation criteria, such as EUCAST Expert Rules (3) or CLSI M07 (4), should be used. It is necessary to take following parameters into consideration when interpreting results: species identification, sample origin, patient case history, or results of additional tests.
Another option of evaluation is to use the expert system in the ErbaExpert program.

Quality control: We recommended all following control strains for internal testing of functionality of the antibiotics in the laboratory. Follow EUCAST or CLSI standards when evaluating results. Fresh strains must be used for quality control.

CCM 3955 <i>Pseudomonas aeruginosa</i> (ATCC 27853)											
MIC (mg/l)											
AMS	PIP	PIT	CAZ	AZT	MER	GEN	AMK	COL	CIP	TGC	T/S
-	1–8	1/4–8/4	1–4	2–8	0.125–1	0.5–2	1–4	0.5–2	0.125–1	-	-

CCM 3954 <i>Escherichia coli</i> (ATCC 25922)											
MIC (mg/l)											
AMS	PIP	PIT	CAZ	AZT	MER	GEN	AMK	COL	CIP	TGC	T/S
2/1–8/4	1–4	1/4–8/4	0.06–0.5	0.06–0.5	0.008–0.06	0.25–1	0.5–4	0.25–1	0.004–0.015	0.03–0.25	≤ 0.5/9.5

CCM 4225 <i>Escherichia coli</i> (ATCC 35218)											
MIC (mg/l)											
AMS	PIP	PIT	-	-	-	-	-	-	-	-	-
8/4–32/16	> 64	0.5/4–2/4									

NCTC 13846 <i>Escherichia coli</i> (mcr-1 positive)											
MIC (mg/l)											
-	-	-	-	-	-	-	COL 4* (exceptionally 2–8)	-	-	-	-

CCM 4987 <i>Klebsiella pneumoniae</i> (ATCC 700603)											
MIC (mg/l)											
-	PIT 8–32	-	-	-	-	-	-	-	-	-	-

*The value is valid for ≥ 90% isolates and should only on occasion be 2 or 8 mg/L (source: EUCAST)
ATCC – American Type Culture Collection
CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ, Tel. +420 549 491 430, Fax +420 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz
Health protection: Components of the kit are not classified as dangerous.
Disposal of the used material: Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration. Put paper packaging waste to recycling.
Literature:
(1) The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MIC and zone diameters, <http://www.eucast.org>
(2) CLSI: Performance Standards for Antimicrobial Susceptibility Testing; CLSI dokument M100. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
(3) EUCAST expert rules in antimicrobial susceptibility testing, <http://www.eucast.org>
(4) CLSI: Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; CLSI dokument M07. Wayne, PA: Clinical and Laboratory Standards Institute.

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Литература:
(1) The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MIC and zone diameters, <http://www.eucast.org>
(2) CLSI: Performance Standards for Antimicrobial Susceptibility Testing; CLSI документ M100. Wayne, PA: Clinical and Laboratory Standards Institute.
(3) EUCAST expert rules in antimicrobial susceptibility testing, <http://www.eucast.org>
(4) CLSI: Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; CLSI документ M07. Wayne, PA: Clinical and Laboratory Standards Institute.

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