## Suspensní médium MIC

Kat. č.: MLT00070

Pro mikrobiologii

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Cat. No.: MLT00070

For microbiology

Suspension Medium MIC is supplementary preparation for MIC kits (MIC G-I, G-II, URINE, NEFERM, STAPHY), which are designed for antibiotic susceptibility testing.

The Suspension Medium MIC is supplied in glass tubes; the content of one tube (13 ml) is sufficient for preparation of the suspension of one strain (inoculation of one MIC plate).

The preparation of the suspension should be carried out according to the procedure recommended for a particular type of MIC test.

### Content:

30 tubes with suspension medium MIC  
Instructions for use

### Storage, expiration:

the Suspension medium for MIC should be stored at the temperature of (2–25) °C. Store away from direct sunlight. Date of expiry is indicated on each package.

### Quality control:

Prepare suspension of control strain with density of McF 0.5. Transfer 60 µl of suspension into 13 ml of suspension medium MIC in tube. Incubate at 35±2 °C for 16–20 hours.

*Escherichia coli* CCM 3954 (ATCC 25922) – turbid growth

*Pseudomonas aeruginosa* CCM 3955 (ATCC 27853) – turbid growth

*Staphylococcus aureus* CCM 4223 (ATCC 29213) – turbid growth

### Warning and precautions:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

### Waste disposal:

The residue of the suspension medium should be disposed in accordance with local and national regulations related to the safe handling of such materials. Recycle packaging paper waste.

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