BACITRACIN S

Kat. č.: MLT00084

Pro mikrobiologii

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Cat. No.: MLT00084

For microbiology

Diagnostic discs BACITRACIN S are intended for simple routine presumptive determination of β -hemolytic group A streptococci and for distinguishing of staphylococci and macrococci from other catalase positive cocci.

Principle:

β -hemolytic group A streptococci and some catalase positive cocci are highly susceptible to low concentration of bacitracin (0,04 IU), thus creating defined inhibition zone around the disc BACITRACIN S. Other β -hemolytic streptococci, staphylococci and macrococci are resistant to this concentration or show only very small zone of inhibition.

A package contains:

- 100 diagnostic discs for 100 determinations
- working leaflet

Storage, expiration:

Discs BACITRACIN S should be stored in dark in a refrigerator at (+2 to +8)°C until the expire date indicated on the package.

Required material:

- Blood agar plates
- common laboratory equipment (loops, sterile forceps, burner)

Recommended procedure:

- Streak blood agar plates with loop or swab so as to obtain confluent growth.
- Aseptically place disc BACITRACIN S on the surface of the medium
- Incubate in an atmosphere containing (5–10) % CO₂ at 37°C for (18–24) h.

Caution:

- For professional use only.
- After use all discs must be autoclaved or incinerated.

Evaluation:

After overnight incubation read the diameters of inhibition zones:

Inhibition zone (mm)	Result
> 10	β -hemolytic group A streptococci, micrococci, Kocuria spp.
< 10	β -hemolytic streptococci of other serological group, staphylococci, macrococci

Note:

Some strains of β -hemolytic streptococci of group C and G may show false positive reaction, therefore further characteristics must be taken in consideration (CAMP test, Lancefield grouping etc.).

Standard control strains:

Quality control of diagnostic discs BACITRACIN S may be provided with standard control strains:

Strain	Inhibition zone (mm)
Streptococcus pyogenes CCM 4425	> 10
Kytococcus sedentarius CCM 2699	> 9
Streptococcus CCM 4685	0

These strains are distributed by:

CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ

Tel. +420 549 491 430, Fax +420 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz

Strains are supplied in freeze-dried form or on gelatin discs.

Date of revision: 18. 11. 2020

USED SYMBOLS

REF Catalogue number **IVD** In vitro diagnostics Manufacturer See instruction for use

LOT Lot number Storage temperature Expiry date

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbalachema.com

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