NOVOBIOCIN

Kat. č.: MLT00085

Pro mikrobiologii

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Cat. No.: MLT00085

For microbiology

Diagnostic discs NOVOBIOCIN are intended for differentiation of coagulase negative staphylococci by means of growth test.

Principle:

Test is based on the natural resistance of some species of coagulase negative staphylococci to novobiocin.

A package of NOVOBIOCIN contains:

- 100 diagnostic discs for 100 determinations
- working leaflet

Storage, expiration:

Discs NOVOBIOCIN should be stored in dark in a refrigerator at (+2 to +8)°C until the expire date indicated on the package.

Required material:

- Petri dishes with Mueller-Hinton agar
- common laboratory equipment (loops, sterile forceps, burner)
- McFarland turbidity standards or densitometer

Recommended procedure:

- Prepare 3-hour broth culture of tested strain with turbidity adjusted to value of 0,5 McFarland.
- Overlay Mueller-Hinton agar with inoculum, streak uniformly and pre-dry the plates if necessary.
- Aseptically place disc NOVOBIOCIN.
- Incubate 24 h at 37°C.

Caution:

- For professional use only.

Evaluation:

After incubation measure the diameter of inhibition zones:

Inhibition zone (mm)	Resistance to novobiocin
< 16	resistant
≥ 16	sensitive

Resistant species of coagulase negative staphylococci:

Human – Staphylococcus hominis subsp. novobiosepticus, S. xylosus, S. saprophyticus subsp. saprophyticus, S. cohnii subsp. cohnii, S. cohnii subsp. urealyticum, S. sciuri s. sciuri, S. lentus.

Veterinary – S. gallinarum, S. arlettae, S. equorum, S. kloosii, S. sciuri s. rodentium, S. sciuri s. carnaticus, S. vitulinus a S. saprophyticus s. bovis.

The other up to date described species of coagulase negative staphylococci are sensitive to novobiocin.

Waste disposal:

Used disc should be treated as potentially infectious and should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging waste to recycling or in municipal waste.

Standard control strains:

Quality control of diagnostic discs NOVOBIOCIN may be provided with standard control strains:

Strain	Resistance to novobiocin
Staphylococcus cohnii, ssp. urealyticum CCM 4296	+
Staphylococcus epidermidis CCM 4418	–



These strains are distributed by:

CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ

Tel. +420 549 491 430, Fax +420 549 498 289, http://www.sci.muni.cz/ccm, e-mail: ccm@sci.muni.cz

Strains are supplied in freeze-dried form or on gelatin discs.

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