Kat. č.: MLT00029

mCOLItest

Pro mikrobiologii

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Cat. No.: MLT00029

For microbiology

The mCOLItest is intended for the detection of characters typical for *Escherichia coli* – the activity of β -glucuronidase, for the purpose of microbiological control of samples by applying the membrane filter method. The preparation can be used for the control of drinking water, in the food industry etc. Four incubation hours are sufficient for the mCOLItest evaluation. The preparation enables the β -glucuronidase detection with colonies grown upon agar media that are recommended for the detection of coliform, resp. thermotolerant coliform bacteria, such as mFC agar, MacConkey agar, Endo agar. A single package of mCOLItest enabling 15 determinations.

Principle:

The membrane filter with bacterial colonies shall be carried over upon the mCOLItest disc surface and allowed to incubate for 4 hours. Afterwards, evaluate the reaction of β -glucuronidase. In the capacity of the substrate serving for the detection of the β -glucuronidase activity, there is 4 methylumbelliferyl- β -D-glucuronide there, being hydrolyzed under the presence of β -glucuronidase – the liberated 4 methylumbelliferon giving a blue fluorescence under a UV radiation source.

A package of the mCOLItest contains:

- 5 bags with discs serving for the detection of β -glucuronidase and with a desiccant, each bag containing 3 pcs of discs for 3 determinations
- Working leaflet

Caution: For professional use only.

Storage: mCOLItest should be stored at the temperature of (+2 to +8) °C. Store away from direct sunlight. The expiry date is indicated on each package.

Required material:

- Pipette 0.9 ml, sterile physiological solution
- Empty Petri dishes with a diameter of 6 cm
- Tweezers
- Thermostat 37 °C
- UV lamp, Cat. No. 50001471
- Usual equipment of a microbiological laboratory (burner, markers etc.)

Recommended procedure:

- Insert 1 disc from the mCOLItest into the empty Petri dish.
- Moisten the disc with 0.9 ml of the sterile physiological solution.
- Carry carefully the membrane filter with colonies over onto the disc by using tweezers so that between the disc and the filter, there are no air bubbles left there (it is recommended to lay the filter on the disc gradually, from one filter edge to the other one).
- Cover the Petri dish with a lid, insert into the thermostat and incubate for 4 hours at the temperature of (35–37) °C.
(Note: resp. before the incubation, record the position of colonies of the coliform, resp. thermotolerant coliform bacteria on the filter upon the dish lid using the marker, at the same time, remark their mutual position at the very moment of recording the colonies on one side of the Petri dish as well as on the lid).
- Read off the reaction on the mCOLItest after the incubation period has elapsed: evaluate the fluorescence reaction in the dark, under a source of the UV radiation with a wave length of approx. 360 nm.
- **The blue fluorescence** appears with the respective colony in the case of the positive β -glucuronidase reaction.
- The colonies of coliform, resp. thermotolerant coliform bacteria delivering the positive β -glucuronidase reaction should be considered as presumptive *Escherichia coli*.

Note:

- Reading of the fluorescence reaction should be carried out in the dark: in a dark room, on a well shaded working table, preferentially in a device serving for the fluorescence evaluation (e.g. the Hansen's box).
- The fluorescence reaction can be set off with ammonia vapours, if need be: drip 1 drop of the concentrated ammonia upon the lid of the Petri dish and put the dish with the mCOLItest and filter onto the lid, incubate for approx. 10 minutes at the laboratory temperature. Afterwards, read the fluorescence reaction.
- Insert the used material into the vessel intended for the infectious material and autoclave or destroy it by incineration.
- In the case of not using up the bag contents immediately after opening, seal-in the supporting plates left together with the desiccant back into the bag, resp. store them dry and keep in dark, nevertheless not longer than 14 days after having the bag open.

Standard control strains:

For the quality control of the mCOLItest as well as the interpretation of the positive and negative reaction, the recommended control strains can be applied: using the inoculation loop, apply the growth of the control strain from the agar medium directly to the surface of the mCOLItest disc, moistened with 0.9 ml of the sterile saline. Incubate mCOLItest disc and read the reaction according to the Recommended procedure paragraph.

Control strains:

Strain	Reaction
<i>Escherichia coli</i> CCM 3988	positive, blue fluorescence
<i>Serratia marcescens</i> CCM 303	negative, no fluorescence

These strains are supplied in by the Czech Collection of Microorganisms, Masaryk University, Kamenice 5, 602 00 Brno, Czech Republic
Phone: (+420) 549 491 430, Fax: (+420) 543 247 339, e-mail: ccm@sci.muni.cz.
The strains are delivered in a lyophilized form or on gelatinized discs.

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