HIPPURATEtest

Kat. č.: MLT00036

Pro mikrobiologii

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Cat. No.: MLT00036

For microbiology

Detection strip HIPPURATEtest is intended for the determination of sodium hippurate hydrolysis. It is recommended for the presumptive identification of Group B streptococci, *Gardnerella vaginalis* and *Campylobacter jejuni*. One packing of HIPPURATEtest enables 50 determinations.

Principle:
Sodium hippurate is split by enzyme hippurate hydrolase into glycine and benzoic acid.
Liberated glycine is detected by colour reaction with ninhydrin, which is contained in the Reagent for the test HIPPURATE.

Sodium hippurate + H₂O Hippurate hydrolase Benzoic acid + Glycine

A package of the HIPPURATEtest contains:
• 50 detection strips for 50 determinations
• working instructions

Reagent:
Reagent for HIPPURATEtest (Cat. No.: MLT00019)
1 package of the reagent is sufficient for 4 packages of the HIPPURATEtest.

Storage:
HIPPURATEtest should be stored at the temperature of (+2 to +8) °C. Store away from direct sunlight. The expiry date is indicated on each package.

Recommended procedure:
• Use the pure culture of the tested strain (24 hours old culture)
• Take a tube (100x15 mm) with 0.5 ml of sterile saline Emulsify the pure culture with a loop in the tube by rubbing against the side of the tube and mixing thoroughly with the saline
• The density of the suspension should be equal to number 2 of McFarland's turbidity scale
• Insert the strip into the tube – the zone of the strip has to be soaked into the suspension
• Incubate at 37 °C for 22–24 hours
• Add 4 drops of the Reagent for HIP into the tube after time of incubation
• Incubate 5–10 minutes at the temperature of (+15 to +25) °C
• Read the reaction visually in accordance with the table „Evaluation of the sodium hippurate hydrolysis test“.

Caution:
• For professional use only

Waste disposal:
Used strip should be treated as potentially infectious and should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging waste to recycling or in municipal waste.

Test	Code	Reaction	
Sodium hippurate	HIP	positive	negative
		blue-violet to blue ring	colourless to bluish

It is recommended to use standard control strains for checking of the reaction and interpretation of the colours of the reaction.

Strain	CCM	HIP
<i>Streptococcus uberis</i>	4617	positive
<i>Kytococcus sedentarius</i>	2699	negative

These strains are distributed by
Czech Collection of Microorganisms, Masaryk University,
Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ
Tel.: +420 549 491 430, Fax: +420 549 498 289 <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz
The strains are supplied lyophilised or on gelatin discs.

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