## Suspennzí médium MIC G+ 200

(CZ)

Kat. č.: MLT00075

Pro mikrobiologii

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Cat. No.: MLT00075

For microbiology

Suspension Medium MIC G+ 200 is supplementary preparation for MIC G+ kit which is designed for antibiotic susceptibility testing. Suspension Medium MIC G+ 200 is supplied in glass bottles. The preparation of the suspension should be carried out according to the procedure recommended for MIC G+ test. The preparation of suspension is also included in the instructions for ErbaFill.

### Content:

2 bottles with suspension medium MIC G+ (200 ml)  
Instructions for use

### Storage, expiration:

Suspension medium for MIC G+ should be stored at the temperature of (2–25) °C. Store away from direct sunlight. Date of expiry is indicated on each package.

### Quality control:

Aseptically transfer 5 ml of Suspension medium G+ 200 into tube. Prepare suspension of control strain with density of McF 0.5. Transfer 60 µl of suspension into suspension medium MIC G+ in tube. Incubate at 35±2 °C for 16–20 hours.  
*Enterococcus faecalis* CCM 4224 (ATCC 29212) – turbid growth

### Warning and precautions:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

**Hazards identification in accordance with Regulation (EC) No 1272/2008**  
Reagent is not classified as dangerous.

### Waste disposal:

The residue of the suspension medium should be disposed in accordance with local and national regulations related to the safe handling of such materials. Recycle packaging paper waste.

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