V+K DISK

Kat. č.: MLT00087

Pro mikrobiologii

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Cat. No.: MLT00087

For microbiology

Diagnostic discs V+K DISK are intended for the selective isolation of *Neisseria meningitidis*.**Principle:**

In the place of diffusion of antibiotics (vancomycin+colistin) from disc the zone of partially selective medium is created, thus inhibiting the growth of accompanying flora.

A package of V+K DISK contains:

- 100 diagnostic discs for 100 determinations
- working leaflet

Supplementary material available from Erba Lachema s.r.o.:

- NEISSERIAtest, Cat. No. MLT00011 – for 36 determinations
- Suspension medium for NEISSERIAtest, Cat. No. MLT00025 – for 18 determinations
- ONPtest, Cat. No. MLT00038 – for 50 determinations
- OXitest, Cat. No. MLT00039 – for 50 determinations

Storage, expiration:

Discs V+K DISK should be stored in dark in a refrigerator at (+2 to +8) °C until the expire date indicated on the package.

Required material:

- Kits and material supplied by Erba Lachema s.r.o.
- Blood agar plates
- common laboratory equipment (loops, sterile forceps, tubes, burner)

Recommended procedure:

- Streak blood agar plate heavily with nasopharyngeal swab or with other clinical material.
- Aseptically place disc V+K DISK on the surface of the medium
- Incubate in an atmosphere containing (5–10) % CO₂ at 37 °C for (18–24) h.

Note:

Suitable combination of diagnostic disks (V+K DISK, BACITRACIN 10 UI, OPTOCHIN) may provide the determination of majority of important pathogens from upper respiratory tract in a single sample.

Evaluation:In the zone of diffusion of antibiotics grow *Neisseria meningitidis*, *Neisseria lactamica* and yeasts. Their differentiation:

- yeasts – by means of direct microscopic examination in Gram-stained smear
- *Neisseria lactamica* – by means of ONPtest – positive result up to 1 h
- *Neisseria meningitidis* – identification by means of NEISSERIAtest, result up to 4 h – possibility of serotyping of the strain at the same day

Cautions:

- For professional use only.
- After use all discs must be autoclaved or incinerated.

Standard control strains:

Quality control of diagnostic discs V+K DISK may be provided with standard control strains:

Strain	Growth in the zone of antibiotics
<i>Neisseria meningitidis</i> CCM 4634	growth
<i>Streptococcus pneumoniae</i> CCM 4501	inhibition of the growth > 20 mm

These strains are distributed by:

CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science,
Kamenice 5, building A25, 625 00 Brno, CZ
Tel. +420 549 491 430, Fax +420 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz
Strains are supplied in freeze-dried form or on gelatin discs.

Date of revision: 18. 11. 2020

USED SYMBOLS

REF	Catalogue number	IVD	In vitro diagnostics	Manufacturer	See instruction for use
LOT	Lot number	Storage temperature	Expiry date		

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbalachema.com

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