

ADVISORY NOTICE

LAURA *version 2* Last Result presentation

LAURA *version 2*, Catalogue No.: INS00063


Dear customer,

we would like to inform you that the latest SW version (1.13) of LAURA *version 2* includes the possibility, which allows the last measured result to be displayed on the main screen during the measurement. Enabling this function could be switched on by user in the customization submenu (Menu → Setup → Customization → Last result).



In case of any questions, please contact your local ERBA team.

Best regards,



Erba Lachema s.r.o.

Technical Product Management

Brno, 01th February 2022

Erba Lachema s.r.o.

Karásek 1d, 621 33 Brno

DIČ: CZ26918846 -2-