

ADVISORY NOTICE

EC cartridge Biosensors

Catalogue no.	Product name
REG00056	EC cartridge M
REG00057	EC cartridge L
REG00058	EC cartridge S
REG00061	EC cartridge S plus iCa
REG00062	EC cartridge M plus iCa
REG00063	EC cartridge L plus iCa
REG00064	EC Biosensor
REG00065	EC cartridge XS
REG00066	EC cartridge XS plus iCa
REG00067	EC cartridge S/60D
REG00068	EC cartridge S plus iCa/60D
REG00069	EC cartridge XS/60D
REG00070	EC cartridge XS plus iCa/60D

Dear valued customer,

to help provide better environmental control of *EC cartridge Biosensors* during transportation, we would like to inform you that starting 27th June 2022, all biosensors will be individually foil-packed and included in the cartridge pack.

To prepare the cartridge for use, all that is required is to remove the Biosensor from the foil and insert it into EC Cartridge as described in the user manual.

The package consisting of EC cartridge and separate Biosensor must be shipped at chilled transport conditions (2 - 8°C).

We thank you in advance for your understanding – in case of any questions, please, contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.

Technical product management

Brno 17th August 2022


Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Ref.No.AN_20_2022