

ASSAY VALUES AND EXPECTED RANGES

Atestované hodnoty a očekávaná rozmezí/Valores de la media y rangos esperados

QCP DATA MONTHS/Použitelnost/ Fecha de la caducidad:
MONTH IN A YEAR/měsíce roku/ mes del año
JULY, AUGUST, SEPTEMBER / červenec, srpen, září
LOT
U0717

05-09-2017

Name/Název/ Nombre	Cat.No./Kat.č./No.Cat.	Level/Hladina/ Nivel	Package volume/ Objem balení/Volumen
ELite 5 HEM Control Low	HEM00009	Low	3 ml
ELite 5 HEM Control Normal	HEM00010	Normal	3 ml
ELite 5 HEM Control High	HEM00011	High	3 ml


IVD


Before using refer to the instruction sheet for mixing directions. Před použitím čtěte návod. Lea las instrucciones antes de usar.

Instruments/Přístroj/Instrumento		ERBA Elite 5		
Parameter/ Analyt/ Analito		CONTROL "Low"	CONTROL "Normal"	CONTROL "High"
		LOT:U0717L	LOT:U0717N	LOT:U0717H
		Mean ± Limit	Mean ± Limit	Mean ± Limit
WBC	x 10 ³ /μL	3.40±0.50	8.60±1.00	23.20±2.50
NEU	x 10 ³ / μL	1.33±0.66	4.75±1.89	14.60±5.20
LYM	x 10 ³ / μL	1.62±0.70	2.80±1.48	5.15±3.38
MONO	x 10 ³ / μL	0.39±0.39	0.87±0.87	2.92±2.92
EOS	x 10 ³ / μL	0.03±0.03	0.09±0.09	0.30±0.30
BASO	x 10 ³ / μL	0.03±0.03	0.09±0.09	0.23±0.23
NEU	%	39.0±12.0	55.2±14.0	62.9±14.0
LYM	%	47.6±12.0	32.6±12.0	22.2±11.0
MONO	%	11.4±11.4	10.1±10.1	12.6±12.6
EOS	%	1.0±1.0	1.1±1.1	1.3±1.3
BASO	%	1.0±1.0	1.0±1.0	1.0±1.0
RBC	x 10 ⁶ / μL	2.20±0.18	4.66±0.24	5.14±0.30
HGB	g/dL	5.9±0.4	13.1±0.5	15.8±0.7
HCT	%	18.6±2.0	43.6±2.5	51.7±3.0
MCV	fL	84.5±5.0	93.5±5.0	100.5±5.0
MCH	pg	26.8±2.8	28.1±2.8	30.7±2.8
MCHC	g/dL	31.7±3.6	30.1±3.0	30.6±3.0
RDW-SD	fL	50.0±10.0	51.0±10.0	52.0±10.0
RDW-CV	%	17.6±5.0	15.8±5.0	14.0±5.0
PLT	x 10 ³ / μL	70±25	240±45	450±70
PDW-CV	%	42.0±9.9	41.5±9.9	40.0±9.9
PDW-SD	fL	20.0±5.0	19.0±5.0	19.5±5.0
MPV	fL	7.5±3.0	8.0±3.0	8.3±3.0
PCT	%	0.05±0.04	0.19±0.08	0.38±0.10
Gravity X			125	
Gravity Y			81	

