


DENSILAMETER® II



**Comfortable Determination
of Microbial Suspension Density**



DENSILAMETER® II

Patented Method of Inoculums Measurement

- Elimination of subjective errors related to determination of inoculum density
- Economic efficiency
- Measurement accuracy
- Measurement of turbidity from 0.0 to 15.0 McFarland
- Standard and user calibration available
 - user calibration as optional
- Variability of user test tubes
- Low volume of suspension needed



Specifications are subject to change without prior notice.
Pictures used only for representation purposes - subject to change without prior notice.

Erba Lachema s.r.o.
 Karásek 2219/1d, 621 00 Brno, Czech Republic
 Tel.: +420 517 077 111
 E-mail: diagnostics@erbamannheim.com
www.erbalachema.com



Contact info / map

DLMII20180901EN