

MAGO 4

Allergy & food intolerance



- Automated processing of ELISA, IFA and agglutination
- MAGO 4 pre-programmed test procedures
- Up to 12 ELISA or 16 IFA simultaneously in one batch
- Wide range of ELISA, IFA and ImmunoBlot tests

Allergy & food intolerance

Allergy

ALKS-100	Specific IgE Base kit	96 tests
	Specific IgE - Basic kit	
ALKC-102	IgE Standard Curve	8 standards
	Specific IgE - Standard with curves	
ALKT-103	Total IgE	10 x 6 tests
	Total IgE kit	
AL	Allergens	3 ml / 24 tests
	280 liquid allergens and 17 mixtures	

Food intolerance

G-090	IgG Screen Nutritional kit	3 x 96 tests
	ELISA test for semi-quantitative analysis of antibodies to 90 food allergens in human serum	
DBL-01	β-Lactoglobulin IgG	96 tests
	ELISA test for quantitative detection of IgG antibodies to cow's milk β -Lactoglobulin in human serum	
DBL-10	Soya IgG	96 tests
	ELISA test for quantitative detection of IgG antibodies to Soya protein	

Celiac disease

I-528	Anti-Tissue Transglutaminase IgA	96 tests
	Test with human recombinant antigen for the IgA anti-human tTG	
I-527	Anti-Tissue Transglutaminase IgG	96 tests
	Test with human recombinant antigen for the IgG anti-human tTG	
I-526	Anti-Gliadin IgA	96 tests
	Test for the IgA anti-Deamidated Gliadin	
I-524	Anti-Gliadin IgG	96 tests
	Test for the IgG anti-Deamidated Gliadin	
I-522	Anti-Parietal Cell	96 tests
	Test for the anti-Parietal Cells Abs	
I-513	Anti-Intrinsic Factor Antibody ELISA	96 tests
	Test for the anti-Intrinsic Factor Abs	
F-127A	CELIAC - Primate Distal Esophagus (EMA) IFA Kit	8 x 6 tests
	EMA test for the IgA with primate distal esophagus	

Complete list of products in Delta Biologicals catalogue.

