

SYSPACKS™

Reagents for ERBA XL Analysers



 **Erba Lachema s.r.o.**
Karásek 2219/1d, 621 00 Brno, Czech Republic
Tel.: +420 517 077 111
E-mail: diagnostics@erbamannheim.com
www.erbalachema.com



Contact info / map

Ready-to-Use Reliable

BLT20180901EN

Specifications are subject to change without prior notice.
Pictures used only for representation purposes - subject to change without prior notice.

CLINICAL CHEMISTRY

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0001	ALBUMIN	ALB 440	R1:10 x 44 ml			BCG
XSYS0002	ALKALINE PHOSPHATASE	ALP 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. IFCC (AMP buffer)
XSYS0003	ALPHA AMYLASE (SINGLE REAGENT)	AMY 110	R1: 5 x 22 ml			CNP-G, single reagent
XSYS0017	ALT/GPT	ALT/GPT 330	R1: 6 x 44 ml	R2: 3 x 22 ml		mod. IFCC
XSYS0016	AST/GOT	AST/GOT 330	R1: 6 x 44 ml	R2: 3 x 22 ml		mod. IFCC
XSYS0028	BILIRUBIN DIRECT	BIL D 330	R1: 6 x 44 ml	R2: 3 x 22 ml		Walter-Gerarde
XSYS0023	BILIRUBIN TOTAL	BIL T 330	R1: 6 x 44 ml	R2: 3 x 22 ml		Walter-Gerarde
XSYS0086	BILIRUBIN DIRECT DCA	BIL D DCA 330	R1: 6 x 44 ml	R2: 3 x 22 ml		with DCA
XSYS0087	BILIRUBIN TOTAL DCA	BIL T DCA 330	R1: 6 x 44 ml	R2: 3 x 22 ml		with DCA
XSYS0007	CALCIUM	CA 120	R1: 10 x 12 ml			Arsenaso III
XSYS0022	CREATINE KINASE	CK 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. IFCC,DGKCH
XSYS0029	CREATINE KINASE MB	CK MB 110	R1: 2 x 44 ml	R2: 2 x 11 ml		immunoinhibition
XSYS0024	CREATININE	CREA 275	R1: 5 x 44 ml	R2: 5 x 11 ml		Jaffé
XSYS0085	CREATININE ENZYMATIC	CREA ENZ 200	R1: 5 x 30 ml	R2: 5 x 10 ml		enzymatic
XSYS0011	GAMMAGLUTAMYLTRANSFERASE	GGT 110	R1: 2 x 44 ml	R2: 2 x 11 ml		Persijn van der Slik
XSYS0012	GLUCOSE	GLU 440	R1:10 x 44 ml			GOD,POD
XSYS0043	HDL DIRECT	HDL C 160	R1: 4 x 30 ml	R2: 4 x 10 ml		direct method with PVS,PEGME
XSYS0008	CHLORIDE	CL 120	R1: 10 x 12 ml			thiocyanate
XSYS0009	CHOLESTEROL	CHOL 440	R1:10 x 44 ml			CHOD,PAP
XSYS0049	IRON	IRON 125	R1: 4 x 25 ml	R2: 2 x 12,5 ml	Standard: 2 x 2 ml	ferrozine
XSYS0013	LACTATEDEHYDROGENASE-P	LDH 110	R1: 2 x 44 ml	R2: 2 x 11 ml		optim. DGKCH (pyruvate)
XSYS0044	LDL DIRECT	LDL C 80	R1: 2 x 30 ml	R2: 2 x 10 ml		direct method with PVS,PEGME
XSYS0081	LIPASE	LIP 110	R1: 2 x 44 ml	R2: 2 x 11 ml		enzymatic, colorimetric
XSYS0040	MAGNESIUM	MG 88	R1: 2 x 44 ml			xylidyl blue
XSYS0027	MICROPROTEIN	MP 120	R1: 10 x 12 ml		Standard: 1x 5 ml	pyrogallol red
XSYS0015	PHOSPHORUS	PHOS 120	R1: 10 x 12 ml			UV molybdate
XSYS0018	TOTAL PROTEIN	TP 440	R1:10 x 44 ml			biuret, single reagent
XSYS0041	TRIGLYCERIDES	TG 440	R1:10 x 44 ml			GPO,PAP, single reagent
XSYS0050	UIBC	UIBC 125	R1: 4 x 25 ml	R2: 2 x 12,5 ml	Standard: 2 x 2 ml	ferrozine
XSYS0020	UREA	UREA 275	R1: 5 x 44 ml	R2: 5 x 11 ml		Urease, GLDH
XSYS0021	URIC ACID	UA 275	R1:10 x 44 ml			enzymatic, single reagent
XSYS0042	URIC ACID (SINGLE REAGENT)	UA 440	R1: 5 x 44 ml	R2: 5 x 11 ml		enzymatic, two reagents

IMMUNOTURBIDIMETRIC KITS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE			METHOD
XSYS0046	ASO	ASO	R1: 2 x 40 ml	R2: 2 x 10 ml		latex enhanced immunoturbidimetry
XSYS0047	CRP	CRP	R1: 2 x 40 ml	R2: 2 x 10 ml		immunoturbidimetry
XSYS0048	RF	RF	R1: 2 x 40 ml	R2: 2 x 8 ml		immunoturbidimetry
XSYS0054	HbA1C	HbA1C	R1: 1 x 24 ml	R2a: 1 x 8 ml R2b: 1 x 4 ml	R3: 2 x 50 ml	latex enhanced immunoturbidimetry
XSYS0083	MICROALBUMIN	MAL	R1: 2 x 30 ml	R2: 1 x 10 ml		immunoturbidimetry
XSYS0084	C-REACTIVE PROTEIN HIGH SENSITIVE	CRP-HS	R1: 2 x 40 ml	R2: 2 x 8 ml		immunoturbidimetry

CALIBRATOR, CONTROLS AND OTHERS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE		METHOD
XSYS0034	XL MULTICAL	XL MULTICAL	R1: 4x3 ml		certified human multicalibrator
BLT00080	ERBA NORM	ERBA NORM	R1: 4x5 ml	Water 1x20 ml	certified human multicontrol sera
BLT00081	ERBA PATH	ERBA PATH	R1: 4x5 ml	Water 1x20 ml	certified human multicontrol sera
XSYS0061	HDL/LDL CALIBRATOR	HDL/LDL CAL	R1: 2x1 ml		certified human calibrator for HDL and LDL cholesterol

IMMUNO CALIBRATORS AND CONTROLS FOR SYSTEM PACKS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE		METHOD
XSYS0051	ASO CALIBRATOR	ASO CAL SH	R1: 1 x 1 ml		certified human calibrator for ASO
XSYS0052	RF CALIBRATOR	RF CAL SH	R1: 1 x 1 ml		certified human calibrator for RF
XSYS0053	CRP CALIBRATOR	CRP CAL SH	R1: 1 x 1 ml		certified human calibrator for CRP
XSYS0088	CRP-HS CALIBRATOR	CRP-HS CAL	R1: 1 x 1 ml		certified human calibrator for CRP-HS
XSYS0057	HbA1c CALIBRATOR SET	HBA1C CAL SET	R1: 4 x 0,25 ml		certified human calibrator for HbA1C
BLT20032	MAL CALIBRATOR	MAL CAL	R1: 1 x 1 ml		certified human calibrator for MAL
XSYS0055	HbA1c CONTROL	HBA1C CON L	R1: 1 x 0,25 ml		certified human control serum for HbA1C
XSYS0056	HbA1c CONTROL	HBA1C CON H	R1: 1 x 0,25 ml		certified human control serum for HbA1C
BLT20013	CRP CONTROL	CRP CON H	R1: 1 x 1 ml		certified human control serum for CRP
BLT20014	CRP CONTROL	CRP CON L	R1: 1 x 1 ml		certified human control serum for CRP
BLT20039	RF CONTROL	RF CON	R1: 1 x 1 ml		certified human control serum for RF
BLT20033	MAL CONTROL	MAL CON	R1: 1 x 1 ml		certified human control sera for MAL
BLT20004	ASO CONTROL	ASO CON	R1: 1 x 1 ml		certified human control serum for ASO
BLT20034	MULTICONTROL LEVEL 1	MULTICON L1	R1: 1 x 1 ml		certified human multicontrol sera
BLT20035	MULTICONTROL LEVEL 2	MULTICON L2	R1: 1 x 1 ml		certified human multicontrol sera

SOLUTIONS FOR ISE

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE		METHOD
XSYS0068	XL ISE REAGENT PACK	XL ISE REAGENT PACK	R1: 1x520 ml	R2: 1x190 ml	for the calibration (NA+,K+,CL-,LI+)
XSYS0063	XL ISE 4 CLEANING	XL ISE 4 CLEAN	R1: 1x90 ml	R2: 6x0.3 g	pepsin based cleaner
XSYS0067	XL URINE DIL	XL Urine Diluent	R1: 1x52 ml		urine diluent

SOLUTIONS FOR ISE - ONLY 3CHANNELS SYSTEM

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE		METHOD
XSYS0064	XL Calibrant A	XL CAL A	R1: 3x500 ml		for the calibration (NA+,K+,CL-)
XSYS0065	XL Calibrant B	XL CAL B	R1: 1x52 ml		for the calibration (NA+,K+,CL-)
XSYS0062	XL ISE CLEANING	XL ISE CLEAN	R1: 1x52 ml		kinase based cleaner

WASHING SOLUTIONS

CAT.NO.	PRODUCT NAME (ANALYTE)	PACK NAME	PACKAGE SIZE		METHOD
XSYS0066	XL WASH	XL WASH	R1: 4x100 ml		wash solution for daily washing
XSYS0082	ERBA XL AutoWash AC/AL	XL AUTOWASH AC/AL	R1: 5x44 ml	R1:5x44 ml	acid and alkali solution for ocassional washing