

AST-GOT_R2

Creation date	27th May 2015	Version	4.0
Revision date	29th November 2023		

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- 1.1. Product identifier** AST-GOT_R2
Substance / mixture mixture
Number BLT00050, BLT00051, XSYS0016, XSYS0073
Other mixture names
AST/GOT 250, AST/GOT 500, AST/GOT 330, AST/GOT 564 XL-1000
- 1.2. Relevant identified uses of the substance or mixture and uses advised against**
Mixture's intended use
Diagnostic reagent for quantitative in vitro determination of AST/GOT (Aspartate Aminotransferase) in human serum and plasma.
Main intended use
PC-MED-OTH Other medical devices
Secondary uses
PC-TEC-19 Reagents and laboratory chemicals
The use descriptors
PC 21 Laboratory chemicals
Mixture uses advised against
not available
- 1.3. Details of the supplier of the safety data sheet**
Manufacturer
Name or trade name Erba Lachema s.r.o.
Address Karásek 2219/1d , Brno, 62100
Czech Republic
Identification number (CRN) 26918846
VAT Reg No CZ26918846
Phone +420 517 077 111
E-mail msds@erba.com
Web address www.erbalachema.com
Competent person responsible for the safety data sheet
Name Erba Lachema s.r.o.
E-mail msds@erba.com
- 1.4. Emergency telephone number**
European emergency number: 112 112

SECTION 2: Hazards identification

- 2.1. Classification of the substance or mixture**
Classification of the mixture in accordance with Regulation (EC) No 1272/2008
The mixture is not classified as dangerous according to Regulation (EC) No 1272/2008.
Full text of all classifications and hazard statements is given in the section 16.
Most serious adverse physico-chemical effects
Unknown
Most serious adverse effects on human health and the environment
Unknown.
- 2.2. Label elements**
none
- 2.3. Other hazards**
The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605. Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended.

AST-GOT_R2

Creation date 27th May 2015
 Revision date 29th November 2023 Version 4.0

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Mixture of substances specified below.

Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note
Index: 011-004-00-7 CAS: 26628-22-8 EC: 247-852-1	sodium azide	<0,1	Acute Tox. 2, H300 Aquatic Acute 1, H400 (M=1) Aquatic Chronic 1, H410 (M=1) EUH032	1
CAS: 1403-66-3	Gentamicin	0,01	Repr. 1A, H360 STOT RE 1, H372 Aquatic Acute 1, H400 Aquatic Chronic 1, H410	
Index: 613-344-00-7 CAS: 3811-73-2 EC: 223-296-5	Pyridine-2-thiol 1-oxide, sodium salt	0,002	Acute Tox. 4, H302 Acute Tox. 3, H311+H331 Skin Irrit. 2, H315 Skin Sens. 1, H317 Eye Irrit. 2, H319 STOT RE 1, H372 (nervous system) Aquatic Acute 1, H400 (M=100) Aquatic Chronic 2, H411 EUH070 Specific concentration limit: ATE Inhalation (dust/mist) = 0,5 mg/l ATE Dermal = 790 mg/kg bw ATE Oral = 500 mg/kg bw	

Notes

1 A substance for which exposure limits are set.

Full text of all classifications and hazard statements is given in the section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures**

There are not any particular first aid measures required.

If inhaled

Discontinue the exposure, remove casualty to fresh air.

If on skin

After contact with skin, wash immediately with soap and water.

If in eyes

Rinse an open eye (hold eyelids with fingers) with plenty of water for about 15 minutes.

If swallowed

Rinse mouth with water, drink 1/2 l of lukewarm water.

4.2. Most important symptoms and effects, both acute and delayed**If inhaled**

Not expected.

If on skin

Not expected.

If in eyes

Not expected.

If swallowed

Unknown.

AST-GOT_R2

Creation date 27th May 2015
Revision date 29th November 2023 Version 4.0

4.3. Indication of any immediate medical attention and special treatment needed

Symptomatic treatment.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

The mixture is not flammable. The extinguishing agents should be adapted to burning substances in the surrounding area.

Unsuitable extinguishing media

Unknown.

5.2. Special hazards arising from the substance or mixture

None.

5.3. Advice for firefighters

The mixture is not flammable. The measure should be adapted to burning substances in the surrounding area.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment, see Section 8. Observe the principles of work safety in chemical laboratories. Do not eat, drink or smoke.

6.2. Environmental precautions

Due to amount of chemical substances in a mixture, an impact on the environment is not expected.

6.3. Methods and material for containment and cleaning up

Absorb spilled agent with a suitable inert material (sand, earth, vapex) and store contaminated material in containers for collection of hazardous waste. For waste disposal, see Section 13.

6.4. Reference to other sections

See the Section 7, 8 and 13.

SECTION 7: Handling and storage**7.1. Precautions for safe handling**

Observe the principles of safety work in chemical laboratories. Observe the normal operating procedures for handling chemical substances and mixtures. Do not eat, drink or smoke. Use personal protective equipment, see Section 8.

7.2. Conditions for safe storage, including any incompatibilities

Store in tightly closed containers in a cool, dry place intended for this purpose.

Storage temperature min 2 °C, max 8 °C

7.3. Specific end use(s)

Souprava je určena pro in vitro diagnostické prostředky.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters**

The mixture contains substances for which occupational exposure limits are set.

Czech Republic**Government Regulation 195/2021 Coll.**

Substance name (component)	Type	Value	Conversion for ppm	Note
sodium azide (CAS: 26628-22-8)	PEL	0,1 mg/m ³	0,370	skin penetration is significantly involved during exposure, irritating to mucous membranes (eyes, respiratory system) and skin

AST-GOT_R2

Creation date 27th May 2015
Revision date 29th November 2023 Version 4.0

Czech Republic

Government Regulation 195/2021 Coll.

Substance name (component)	Type	Value	Conversion for ppm	Note
sodium azide (CAS: 26628-22-8)	NPK-P	0,3 mg/m ³	0,370	skin penetration is significantly involved during exposure, irritating to mucous membranes (eyes, respiratory system) and skin

European Union

Commission Directive 2000/39/EC

Substance name (component)	Type	Value	Note
sodium azide (CAS: 26628-22-8)	OEL 8 hours	0,1 mg/m ³	Skin
	OEL 15 minutes	0,3 mg/m ³	

8.2. Exposure controls

Not required.

Eye/face protection

Not required.

Skin protection

Not required.

Respiratory protection

Not required.

Thermal hazard

Not known.

Environmental exposure controls

To eliminate the emergency conditions, have pre-prepared a decontamination mixture and appropriate collection vessels for reaction residues. Dispose of reaction residues and decontaminated mixtures as hazardous waste water in accordance with relevant legal regulations.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	liquid
Colour	colourless
Odour	without fragrance
Melting point/freezing point	data not available
Boiling point or initial boiling point and boiling range	data not available
Flammability	data not available
Lower and upper explosion limit	data not available
Flash point	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
pH	9.5-9.7 (undiluted at 20 °C)
Kinematic viscosity	data not available
Solubility in water	data not available
Partition coefficient n-octanol/water (log value)	data not available
Vapour pressure	data not available
Density and/or relative density	data not available
Relative vapour density	data not available
Particle characteristics	data not available
Form	Clear colourless liquid

9.2. Other information

AST-GOT_R2

Creation date	27th May 2015	Version	4.0
Revision date	29th November 2023		

Oxidising properties
Explosive properties

The product has no oxidizing properties.
The product does not have explosive properties.

SECTION 10: Stability and reactivity

10.1. Reactivity

Under normal conditions of use and storage the mixture is stable.

10.2. Chemical stability

The mixture is stable at normal temperature and pressure.

10.3. Possibility of hazardous reactions

Unknown.

10.4. Conditions to avoid

Avoid exposure to heat and solar radiation.

10.5. Incompatible materials

Kyseliny.

10.6. Hazardous decomposition products

Při kontaktu s kyselinami může uvolňovat toxický plyn.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Inhalation of solvent vapors above values exceeding exposure limits for working environment may result in acute inhalation poisoning, depending on the level of concentration and exposure time. No toxicological data is available for the mixture.

Acute toxicity

Based on the available data, the criteria for classification of the mixture are not met.

AST-GOT_R2						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	ATE	28390 mg/kg				Calculation of value
Dermal	ATE	39500000 mg/kg				Calculation of value
Inhalation (vapor)	ATE	150000 mg/l				Calculation of value

Pyridine-2-thiol 1-oxide, sodium salt						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Inhalation (dust/mist)	ATE	0.5 mg/l				
Dermal	ATE	790 mg/kg bw				
Oral	ATE	500 mg/kg bw				

sodium azide						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	27 mg/kg bw				
Inhalation	LC ₅₀	54 mg/m ³	4 hours	Rat		

Skin corrosion/irritation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

AST-GOT_R2

Creation date	27th May 2015		
Revision date	29th November 2023	Version	4.0

Serious eye damage/irritation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Respiratory or skin sensitisation

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Germ cell mutagenicity

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Carcinogenicity

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Reproductive toxicity

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Toxicity for specific target organ - single exposure

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Toxicity for specific target organ - repeated exposure

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

Aspiration hazard

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

11.2. Information on other hazards

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

SECTION 12: Ecological information**12.1. Toxicity**

Data for the mixture are not available. Based on the available data, the criteria for classification of the mixture are not met.

Acute toxicity

sodium azide				
Parameter	Value	Exposure time	Species	Environment
LC ₅₀	680 µg/l		Fish	Fresh water
EC ₅₀ /LC ₅₀	400 µg/l		Invertebrates	Fresh water
EC ₅₀ /LC ₅₀	150 µg/l		Invertebrates	Salt water
EC ₅₀ /LC ₅₀	348 µg/l		Algae	Fresh water
EC ₅₀ /LC ₅₀	5.6 mg/l		Microorganisms	
NOEC	30 µg/l		Microorganisms	

12.2. Persistence and degradability

No data are available for either the mixture or the components.

12.3. Bioaccumulative potential

No data are available for either the mixture or the components.

AST-GOT_R2

Creation date	27th May 2015		
Revision date	29th November 2023	Version	4.0

12.4. Mobility in soil

No data are available for either the mixture or the components.

12.5. Results of PBT and vPvB assessment

Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

12.6. Endocrine disrupting properties

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

12.7. Other adverse effects

Not available.

SECTION 13: Disposal considerations**13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Proceed in accordance with valid regulations on waste disposal. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Empty containers may be used at waste incinerators to produce energy or deposited in a dump with appropriate classification. Perfectly cleaned containers can be submitted for recycling.

Waste management legislation

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

SECTION 14: Transport information**14.1. UN number or ID number**

not subject to transport regulations

14.2. UN proper shipping name

not relevant

14.3. Transport hazard class(es)

not relevant

14.4. Packing group

not relevant

14.5. Environmental hazards

The mixture is not hazardous to the environment during transport.

14.6. Special precautions for user

Reference in the Sections 4 to 8.

14.7. Maritime transport in bulk according to IMO instruments

Not transported.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. Commission Regulation (EU) 2020/878 of 18 June 2020 amending Annex II to Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

15.2. Chemical safety assessment

not available

AST-GOT_R2

Creation date	27th May 2015		
Revision date	29th November 2023	Version	4.0

SECTION 16: Other information

A list of standard risk phrases used in the safety data sheet

H300	Fatal if swallowed.
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H360	May damage fertility or the unborn child.
H372	Causes damage to organs through prolonged or repeated exposure.
H372	Causes damage to the nervous system through prolonged or repeated exposure.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H411	Toxic to aquatic life with long lasting effects.
H311+H331	Toxic in contact with skin or if inhaled.

A list of additional standard phrases used in the safety data sheet

EUH032	Contact with acids liberates very toxic gas.
EUH070	Toxic by eye contact.

Other important information about human health protection

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

Key to abbreviations and acronyms used in the safety data sheet

ADR	European agreement concerning the international carriage of dangerous goods by road
BCF	Bioconcentration Factor
CAS	Chemical Abstracts Service
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures
EC	Identification code for each substance listed in EINECS
EINECS	European Inventory of Existing Commercial Chemical Substances
EmS	Emergency plan
EU	European Union
EuPCS	European Product Categorisation System
IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry
LC ₅₀	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD ₅₀	Lethal dose of a substance in which it can be expected death of 50% of the population
log K _{ow}	Octanol-water partition coefficient
NOEC	No observed effect concentration
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, Bioaccumulative and Toxic
PEL	Permissible Exposure Limit
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Agreement on the transport of dangerous goods by rail

AST-GOT_R2

Creation date	27th May 2015	Version	4.0
Revision date	29th November 2023		

UN	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very Persistent and very Bioaccumulative

Acute Tox.	Acute toxicity
Aquatic Acute	Hazardous to the aquatic environment
Aquatic Chronic	Hazardous to the aquatic environment (chronic)
Eye Irrit.	Eye irritation
Repr.	Reproductive toxicity
Skin Irrit.	Skin irritation
Skin Sens.	Skin sensitization
STOT RE	Specific target organ toxicity - repeated exposure

Training guidelines

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

Recommended restrictions of use

not available

Information about data sources used to compile the Safety Data Sheet

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended.
REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

The changes (which information has been added, deleted or modified)

The version 4.0 replaces the SDS version from 12 April 2021. Changes were made in sections 2, 11, 12, 15 and 16.

Statement

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.