

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Substance / mixture

Number

UFI

Other mixture names

AMY 110, AMY SINGLE 100

ALPHA AMYLASE (SINGLE REAGENT)

mixture

XSYS0003, BLT00006

6RYE-AE3K-3H7A-Q7QR

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Mixture's intended use**

Diagnostic reagent for quantitative in vitro determination of alpha-Amylase in human serum, plasma and urine.

**Main intended use**

PC-MED-OTH

Other medical devices

**Secondary uses**

PC-TEC-19

Reagents and laboratory chemicals

**The use descriptors**

PC 21

Laboratory chemicals

**Mixture uses advised against**

not available

**1.3. Details of the supplier of the safety data sheet****Manufacturer**

Name or trade name

Address

Identification number (CRN)

VAT Reg No

Phone

E-mail

Web address

Erba Lachema s.r.o.

Karásek 2219/1d , Brno, 62100

Czech Republic

26918846

CZ26918846

+420 517 077 111

msds@erba.com

www.erbalachema.com

**Competent person responsible for the safety data sheet**

Name

E-mail

Erba Lachema s.r.o.

msds@erba.com

**1.4. Emergency telephone number**

European emergency number: 112 112

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Classification of the mixture in accordance with Regulation (EC) No 1272/2008**

The mixture is classified as dangerous.

Eye Dam. 1, H318

Full text of all classifications and hazard statements is given in the section 16.

**Most serious adverse effects on human health and the environment**

Causes serious eye damage.

**2.2. Label elements****Hazard pictogram****Signal word**

Danger

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

**Hazardous substances**

potassium thiocyanate

**Hazard statements**

H318 Causes serious eye damage.

**Precautionary statements**

P280 Wear eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

**2.3. Other hazards**

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605. Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended.

**SECTION 3: Composition/information on ingredients****3.2. Mixtures****Chemical characterization**

Mixture of substances and additives specified below.

**Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment**

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note
CAS: 333-20-0 EC: 206-370-1	potassium thiocyanate	<4,4	Acute Tox. 4, H302+H312+H332 Eye Dam. 1, H318 Aquatic Chronic 3, H412 EUH032	
Index: 011-004-00-7 CAS: 26628-22-8 EC: 247-852-1	sodium azide	<0,1	Acute Tox. 2, H300 Aquatic Acute 1, H400 (M=1) Aquatic Chronic 1, H410 (M=1) EUH032	1
Index: 017-013-00-2 CAS: 10043-52-4 EC: 233-140-8	calcium chloride	<0,06	Eye Irrit. 2, H319	1

**Notes**

1 A substance for which exposure limits are set.

Full text of all classifications and hazard statements is given in the section 16.

**SECTION 4: First aid measures****4.1. Description of first aid measures**

Take care of your own safety. If any health problems are manifested or if in doubt, inform a doctor and show him information from this safety data sheet. If unconscious, put the person in the stabilized (recovery) position on his side with his head slightly bent backwards and make sure that airways are free; never induce vomiting. If the person vomits by himself, make sure that the vomit is not inhaled. In life threatening conditions first of all provide resuscitation of the affected person and ensure medical assistance. Respiratory arrest - provide artificial respiration immediately. Cardiac arrest - provide indirect cardiac massage immediately.

**If inhaled**

Terminate the exposure immediately; move the affected person to fresh air.

**If on skin**

Remove contaminated clothes. Wash the affected area with plenty of water, lukewarm if possible. Soap, soap solution or shampoo should be used if there is no skin injury. Provide medical treatment if skin irritation persists.

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

**If in eyes**

Rinse eyes immediately with a flow of running water, open the eyelids (also using force if needed); remove contact lenses immediately if worn by the affected person. No neutralization should be performed in any case! Rinsing should be continued for 10-30 minutes from the inner to the outer eye corner to make sure that the other eye is not involved. Depending on the situation, call medical rescue service or ensure medical treatment as promptly as possible. Everyone must be referred for treatment even if affected only a little.

**If swallowed**

Rinse out the mouth with clean water. In the event of issues, find medical help.

**4.2. Most important symptoms and effects, both acute and delayed****If inhaled**

Inhaling vapours can cause corrosion of the breathing system.

**If on skin**

Not expected.

**If in eyes**

Causes serious eye damage.

**If swallowed**

Corrosion of the digestion system can occur.

**4.3. Indication of any immediate medical attention and special treatment needed**

Symptomatic treatment.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**

Alcohol-resistant foam, carbon dioxide, powder, water spray jet, water mist.

**Unsuitable extinguishing media**

Water - full jet.

**5.2. Special hazards arising from the substance or mixture**

Inhalation of hazardous degradation (pyrolysis) products may cause serious health damage.

**5.3. Advice for firefighters**

Self-Contained Breathing Apparatus (SCBA) with a chemical protection suit only where personal (close) contact is likely. Use a self-contained breathing apparatus and full-body protective clothing. Do not allow run-off of contaminated fire extinguishing material to enter drains or surface and ground water.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment for work. Follow the instructions in the Sections 7 and 8. Prevent contact with skin and eyes. Observe the principles of work safety in chemical laboratories. Do not eat, drink or smoke.

**6.2. Environmental precautions**

Prevent contamination of the soil and entering surface or ground water. Due to amount of chemical substances in a mixture, an impact on the environment is not expected.

**6.3. Methods and material for containment and cleaning up**

Spilled product should be covered with suitable (non-flammable) absorbing material (sand, diatomaceous earth, earth and other suitable absorption materials); to be contained in well closed containers and removed as per the Section 13. In the event of leakage of the substantial amount of the product, inform fire brigade and other competent bodies. After removal of the product, wash the contaminated site with plenty of water. Do not use solvents.

**6.4. Reference to other sections**

See the Section 7, 8 and 13.

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date

01st September 2015

Revision date

16th October 2023

Version

3.0

**SECTION 7: Handling and storage****7.1. Precautions for safe handling**

Prevent formation of gases and vapours in concentrations exceeding the occupational exposure limits. Prevent contact with skin and eyes. Use personal protective equipment as per Section 8. Observe valid legal regulations on safety and health protection. Observe the principles of safety work in chemical laboratories. Observe the normal operating procedures for handling chemical substances and mixtures. Do not eat, drink or smoke.

**7.2. Conditions for safe storage, including any incompatibilities**

Store in tightly closed containers in cold, dry and well ventilated areas designated for this purpose.

Storage temperature

min 2 °C, max 8 °C

**7.3. Specific end use(s)**

Intended for in vitro diagnostic devices.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters**

The mixture contains substances for which occupational exposure limits are set.

**Czech Republic****Government Regulation 195/2021 Coll.**

Substance name (component)	Type	Value	Conversion for ppm	Note
sodium azide (CAS: 26628-22-8)	PEL	0,1 mg/m <sup>3</sup>	0,370	skin penetration is significantly involved during exposure, irritating to mucous membranes (eyes, respiratory system) and skin
	NPK-P	0,3 mg/m <sup>3</sup>	0,370	
calcium chloride (CAS: 10043-52-4)	PEL	2 mg/m <sup>3</sup>		irritating to mucous membranes (eyes, respiratory system) and skin
	NPK-P	4 mg/m <sup>3</sup>		

**European Union****Commission Directive 2000/39/EC**

Substance name (component)	Type	Value	Note
sodium azide (CAS: 26628-22-8)	OEL 8 hours	0,1 mg/m <sup>3</sup>	Skin
	OEL 15 minutes	0,3 mg/m <sup>3</sup>	

**8.2. Exposure controls**

Follow the usual measures intended for health protection at work and especially for good ventilation. This can be achieved only by local suction or efficient general ventilation. Do not eat, drink and smoke during work. Wash your hands thoroughly with water and soap after work and before breaks for a meal and rest.

**Eye/face protection**

Protective goggles or face shield (based on the nature of the work performed).

**Skin protection**

Hand protection: Protective gloves resistant to the product. When choosing appropriate thickness, material and permeability of the gloves, observe recommendations of their particular manufacturer. Observe other recommendations of the manufacturer. Other protection: protective workwear. Contaminated skin should be washed thoroughly.

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

**Respiratory protection**

If all workplace limits are observed and good ventilation is ensured, no special precautions necessary.

**Thermal hazard**

Not available.

**Environmental exposure controls**

Observe usual measures for protection of the environment, see Section 6.2.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state	liquid
Colour	colourless
Odour	without fragrance
Melting point/freezing point	data not available
Boiling point or initial boiling point and boiling range	data not available
Flammability	non-flammable
Lower and upper explosion limit	data not available
Flash point	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
pH	data not available
Kinematic viscosity	data not available
Solubility in water	soluble
Partition coefficient n-octanol/water (log value)	data not available
Vapour pressure	data not available
Density and/or relative density	data not available
Relative vapour density	data not available
Particle characteristics	data not available
Form	Clear colourless liquid

**9.2. Other information**

Oxidising properties	The product has no oxidizing properties.
Explosive properties	The product does not have explosive properties.
none	

**SECTION 10: Stability and reactivity****10.1. Reactivity**

When used in the standard way, there is not any dangerous reaction with other substances.

**10.2. Chemical stability**

The product is stable under normal conditions.

**10.3. Possibility of hazardous reactions**

Unknown.

**10.4. Conditions to avoid**

The product is stable and no degradation occurs under normal use. Protect against flames, sparks, overheating and against frost. Keep away from heat and frost. Store at a temperature between +2°C and +8 °C.

**10.5. Incompatible materials**

Protect against strong acids, bases and oxidizing agents.

**10.6. Hazardous decomposition products**

Not developed under normal uses. May release toxic fumes on contact with acids.

**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008**

Inhalation of solvent vapors above values exceeding exposure limits for working environment may result in acute inhalation poisoning, depending on the level of concentration and exposure time. No toxicological data is available for the mixture.

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date

01st September 2015

Revision date

16th October 2023

Version

3.0

**Acute toxicity**

Based on the available data, the criteria for classification of the mixture are not met.

**ALPHA AMYLASE (SINGLE REAGENT)**

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	ATE	8279 mg/kg				Calculation of value
Dermal	ATE	25150 mg/kg				Calculation of value
Inhalation (vapor)	ATE	251.5 mg/l				Calculation of value

**calcium chloride**

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Dermal	LD <sub>50</sub>	2630 mg/kg		Rat		

**sodium azide**

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD <sub>50</sub>	27 mg/kg bw				
Inhalation	LC <sub>50</sub>	54 mg/m <sup>3</sup>	4 hours	Rat		

**Skin corrosion/irritation**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Serious eye damage/irritation**

Causes serious eye damage. Data for the components of the mixture are not available.

**Respiratory or skin sensitisation**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Germ cell mutagenicity**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Carcinogenicity**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Reproductive toxicity**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Toxicity for specific target organ - single exposure**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**Toxicity for specific target organ - repeated exposure**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date 01st September 2015  
Revision date 16th October 2023 Version 3.0

**Aspiration hazard**

No data are available for either the mixture or the components. Based on the available data, the criteria for classification of the mixture are not met.

**11.2. Information on other hazards**

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

**SECTION 12: Ecological information****12.1. Toxicity**

Data for the mixture are not available. Based on the available data, the criteria for classification of the mixture are not met.

**Acute toxicity**

calcium chloride				
Parameter	Value	Exposure time	Species	Environment
LC <sub>50</sub>	10650 mg/l	96 hours	Fish (Lepomis macrochirus)	
EC <sub>50</sub>	144 mg/l	48 hours	Daphnia (Daphnia magna)	
IC <sub>50</sub>	3130 mg/l	120 hours	Algae	

sodium azide				
Parameter	Value	Exposure time	Species	Environment
LC <sub>50</sub>	680 µg/l		Fish	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	400 µg/l		Invertebrates	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	150 µg/l		Invertebrates	Salt water
EC <sub>50</sub> /LC <sub>50</sub>	348 µg/l		Algae	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	5.6 mg/l		Microorganisms	
NOEC	30 µg/l		Microorganisms	

**12.2. Persistence and degradability**

No data are available for either the mixture or the components.

**12.3. Bioaccumulative potential**

No data are available for either the mixture or the components.

**12.4. Mobility in soil**

No data are available for either the mixture or the components.

**12.5. Results of PBT and vPvB assessment**

Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

**12.6. Endocrine disrupting properties**

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

**12.7. Other adverse effects**

Not available.

**SECTION 13: Disposal considerations**

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

**13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Proceed in accordance with valid regulations on waste disposal. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Empty containers may be used at waste incinerators to produce energy or deposited in a dump with appropriate classification. Perfectly cleaned containers can be submitted for recycling.

**Waste management legislation**

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

**Waste type code**

16 05 06 laboratory chemicals, consisting of or containing hazardous substances, including mixtures of laboratory chemicals \*

**Packaging waste type code**

15 01 10 packaging containing residues of or contaminated by hazardous substances \*

(\*) - Hazardous waste according to Directive 2008/98/EC on hazardous waste

**SECTION 14: Transport information****14.1. UN number or ID number**

not subject to transport regulations

**14.2. UN proper shipping name**

not relevant

**14.3. Transport hazard class(es)**

not relevant

**14.4. Packing group**

not relevant

**14.5. Environmental hazards**

The mixture is not hazardous to the environment during transport.

**14.6. Special precautions for user**

Reference in the Sections 4 to 8.

**14.7. Maritime transport in bulk according to IMO instruments**

not relevant

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. Commission Regulation (EU) 2020/878 of 18 June 2020 amending Annex II to Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

**15.2. Chemical safety assessment**

not available

**SECTION 16: Other information****A list of standard risk phrases used in the safety data sheet**

H300	Fatal if swallowed.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.



**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015	Version	3.0
Revision date	16th October 2023		

H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.
H302+H312+H332	Harmful if swallowed, in contact with skin or if inhaled.

**Guidelines for safe handling used in the safety data sheet**

P280	Wear eye protection.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.

**A list of additional standard phrases used in the safety data sheet**

EUH032	Contact with acids liberates very toxic gas.
--------	--

**Other important information about human health protection**

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

**Key to abbreviations and acronyms used in the safety data sheet**

ADR	European agreement concerning the international carriage of dangerous goods by road
BCF	Bioconcentration Factor
CAS	Chemical Abstracts Service
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures
EC	Identification code for each substance listed in EINECS
EC <sub>50</sub>	Concentration of a substance when it is affected 50% of the population
EINECS	European Inventory of Existing Commercial Chemical Substances
EmS	Emergency plan
EU	European Union
EuPCS	European Product Categorisation System
IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
IC <sub>50</sub>	Concentration causing 50% blockade
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry
LC <sub>50</sub>	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD <sub>50</sub>	Lethal dose of a substance in which it can be expected death of 50% of the population
log K <sub>ow</sub>	Octanol-water partition coefficient
NOEC	No observed effect concentration
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, Bioaccumulative and Toxic
PEL	Permissible Exposure Limit
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Agreement on the transport of dangerous goods by rail
UN	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very Persistent and very Bioaccumulative

**ALPHA AMYLASE (SINGLE REAGENT)**

Creation date	01st September 2015		
Revision date	16th October 2023	Version	3.0

Acute Tox.	Acute toxicity
Aquatic Acute	Hazardous to the aquatic environment
Aquatic Chronic	Hazardous to the aquatic environment (chronic)
Eye Dam.	Serious eye damage

**Training guidelines**

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

**Recommended restrictions of use**

not available

**Information about data sources used to compile the Safety Data Sheet**

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended.  
REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

**The changes (which information has been added, deleted or modified)**

The version 3.0 replaces the SDS version from 20 March 2021. Changes were made in sections 1, 2, 11, 12, 13, 15 and 16.

**Statement**

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.