

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- 1.1. Product identifier**
Substance / mixture
Number
ERBA LYTE PRO PLUS REAGENT PACK_SOL R
mixture
REG00014
- 1.2. Relevant identified uses of the substance or mixture and uses advised against**
mixture's intended use
ERBA LYTE PRO PLUS REAGENT PACK is intended to calibrate the electrodes, it functions as carrying medium. Solution R is used for calibration of pressure sensor. Use only with ERBA LYTE analyzers.
- The use descriptors**
PC 21 Laboratory chemicals
- Mixture uses advised against not available
- 1.3. Details of the supplier of the safety data sheet**
Manufacturer
Name or trade name
Address
Identification number (CRN)
VAT Reg No
Phone
E-mail
Web address
Erba Lachema s.r.o.
Karásek 2219/1d , Brno, 62100
Czech Republic
26918846
CZ36918846
+420 517 077 111
msds@erbamannheim.com
www.erbalachema.com
- Competent person responsible for the safety data sheet**
Name
E-mail
Erba Lachema s.r.o.
msds@erbamannheim.com
- 1.4. Emergency telephone number**
Poisoning information centre, Na Bojišti 1, Praha, Czech Republic, Tel.: non-stop +420 224 919 293 or +420 224 915 402, Information on health risks only - acute poisoning of humans and animals.

SECTION 2: Hazards identification

- 2.1. Substance or mixture classification**
Classification of the mixture in accordance with Regulation (EC) No 1272/2008
The mixture is classified as dangerous.

Eye Dam. 1, H318

Full text of all classifications and hazard statements is given in the section 16.

Most serious adverse physico-chemical effects

Unknown

Most serious adverse effects on human health and the environment

Causes serious eye damage.

- 2.2. Label elements**
Hazard pictogram

**Signal word**

Danger

Hazardous substances

Lactic acid

Hazard statements

H318 Causes serious eye damage.

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

2.3. Other hazards

Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended.

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Mixture of substances and additives specified below.

Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note.
CAS: 50-21-5 EC: 200-018-0	Lactic acid	8.5-9	Skin Irrit. 2, H315 Eye Dam. 1, H318 Specific concentration limit: Eye Dam. 1, H318: $C \geq 3 \%$ Eye Irrit. 2, H319: $1 \% \leq C < 3 \%$ Skin Irrit. 2, H315: $C \geq 10 \%$	

Full text of all classifications and hazard statements is given in the section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures**

Take care of your own safety. If any health problems are manifested or if in doubt, inform a doctor and show him information from this safety data sheet. If unconscious, put the person in the stabilized (recovery) position on his side with his head slightly bent backwards and make sure that airways are free; never induce vomiting. If the person vomits by himself, make sure that the vomit is not inhaled. In life threatening conditions first of all provide resuscitation of the affected person and ensure medical assistance. Respiratory arrest - provide artificial respiration immediately. Cardiac arrest - provide indirect cardiac massage immediately.

Inhalation

Terminate the exposure immediately; move the affected person to fresh air.

Skin contact

Remove contaminated clothes. Wash the affected area with plenty of water, lukewarm if possible. Soap, soap solution or shampoo should be used if there is no skin injury. Provide medical treatment if skin irritation persists.

Eye contact

Rinse eyes immediately with a flow of running water, open the eyelids (also using force if needed); remove contact lenses immediately if worn by the affected person. No neutralization should be performed in any case! Rinsing should be continued for 10-30 minutes from the inner to the outer eye corner to make sure that the other eye is not involved. Depending on the situation, call medical rescue service or ensure medical treatment as promptly as possible. Everyone must be referred for treatment even if affected only a little.

Ingestion

Rinse mouth with water, drink 1/2 l of lukewarm water, seek medical attention immediately.

4.2. Most important symptoms and effects, both acute and delayed**Inhalation**

Not expected.

Skin contact

Not expected.

Eye contact

Causes serious eye damage.

Ingestion

Corrosion of the digestion system can occur.

4.3. Indication of any immediate medical attention and special treatment needed

Symptomatic treatment.

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

The mixture is not flammable. The measures should be adapted to burning substances in the surrounding area.

Unsuitable extinguishing media

For mixture no unsuitable extinguishing media are known.

5.2. Special hazards arising from the substance or mixture

None.

5.3. Advice for firefighters

No special advice.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment, see Section 8. Observe the principles of work safety in chemical laboratories. Do not eat, drink or smoke. Prevent contact with skin and eyes.

6.2. Environmental precautions

Do not discharge into the drains, surface waters and groundwater.

6.3. Methods and material for containment and cleaning up

Absorb spilled agent with a suitable inert material (sand, earth, vapex) and store contaminated material in containers for collection of hazardous waste. For waste disposal, see Section 13.

6.4. Reference to other sections

See the Section 7, 8 and 13.

SECTION 7: Handling and storage**7.1. Precautions for safe handling**

Prevent contact with skin and eyes. Use personal protective equipment as per Section 8. Observe valid legal regulations on safety and health protection.

7.2. Conditions for safe storage, including any incompatibilities

Store in tightly closed containers in cold, dry and well ventilated areas designated for this purpose.

Storage temperature

min 2 °C, max 25 °C

7.3. Specific end use(s)

For in vitro diagnostic devices.

SECTION 8: Exposure controls/personal protection**8.1. Control parameters**

none

8.2. Exposure controls

Do not eat, drink and smoke during work.

Eye/face protection

Protective goggles.

Skin protection

Protective clothing.

Respiratory protection

Not required with adequate ventilation.

Thermal hazard

Not available.

Environmental exposure controls

Observe usual measures for protection of the environment, see Section 6.2.

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Appearance	clear
Physical state	liquid at 20°C
color	colourless
Odour	odourless
Odour threshold	data not available
pH	6,9 (undiluted)
Melting point/freezing point	0 °C
Initial boiling point and boiling range	100 °C
Flash point	data not available
Evaporation rate	data not available
Flammability (solid, gas)	The product is non-flammable.
Upper/lower flammability or explosive limits	
flammability limits	data not available
explosive limits	data not available
Vapour pressure	data not available
Vapour density	data not available
Relative density	data not available
Solubility(ies)	
solubility in water	data not available
solubility in fats	data not available
Partition coefficient: n-octanol/water	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
Viscosity	data not available
Explosive properties	The product does not have explosive properties.
Oxidising properties	The product has no oxidizing properties.

9.2. Other information

Density	data not available
ignition temperature	data not available

SECTION 10: Stability and reactivity**10.1. Reactivity**

Data not available.

10.2. Chemical stability

The product is stable under normal conditions.

10.3. Possibility of hazardous reactions

Unknown.

10.4. Conditions to avoid

The product is stable and no degradation occurs under normal use.

10.5. Incompatible materials

Protect against strong acids, bases and oxidizing agents.

10.6. Hazardous decomposition products

Not developed under normal uses. Dangerous outcomes such as carbon monoxide and carbon dioxide are formed at high temperature and in fire.

SECTION 11: Toxicological information**11.1. Information on toxicological effects**

No toxicological data is available for the mixture.

Acute toxicity

Based on available data the classification criteria are not met.

Skin corrosion/irritation

May cause slight skin irritation.

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

Serious eye damage/irritation

Causes serious eye damage.

Respiratory or skin sensitisation

Based on available data the classification criteria are not met.

Germ cell mutagenicity

Based on available data the classification criteria are not met.

Carcinogenicity

Based on available data the classification criteria are not met.

Reproductive toxicity

Based on available data the classification criteria are not met.

Toxicity for specific target organ - single exposure

Based on available data the classification criteria are not met.

Toxicity for specific target organ - repeated exposure

Based on available data the classification criteria are not met.

Aspiration hazard

Based on available data the classification criteria are not met.

SECTION 12: Ecological information**12.1. Toxicity****Acute toxicity**

The mixture is not classified as toxic to environment.

12.2. Persistence and degradability

Data not available.

12.3. Bioaccumulative potential

Not available.

12.4. Mobility in soil

Not available.

12.5. Results of PBT and vPvB assessment

Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

12.6. Other adverse effects

Not available.

SECTION 13: Disposal considerations**13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Proceed in accordance with valid regulations on waste disposal. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Empty containers may be used at waste incinerators to produce energy or deposited in a dump with appropriate classification. Perfectly cleaned containers can be submitted for recycling.

Waste management legislation

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

Waste type code

16 05 06 laboratory chemicals, consisting of or containing dangerous substances, including mixtures of laboratory chemicals

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

Packaging waste type code

15 01 10 packaging containing residues of or contaminated by dangerous substances

SECTION 14: Transport information**14.1. UN number**

Not subject to ADR.

14.2. UN proper shipping name

not available

14.3. Transport hazard class(es)

not available

14.4. Packing group

not available

14.5. Environmental hazards

not available

14.6. Special precautions for user

Reference in the Sections 4 to 8.

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not transported.

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. The Act No. 350/2011 Coll., on Chemical Substances and Chemical Preparations as amended (the Chemical Act). Decree No. 432/2003 Coll., laying down conditions for assigning categories to individual jobs, limit values of indices from biological exposure tests, conditions for the sampling of biological materials for biological exposure and the particulars of the reports on work with asbestos and biological agents as amended.

15.2. Chemical safety assessment

Assessment was not carried out.

SECTION 16: Other information**A list of standard risk phrases used in the safety data sheet**

H315 Causes skin irritation.

H318 Causes serious eye damage.

Guidelines for safe handling used in the safety data sheet

P280 Wear protective gloves/protective clothing/eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Other important information about human health protection

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

Key to abbreviations and acronyms used in the safety data sheet

ADR European agreement concerning the international carriage of dangerous goods by road

BCF Bioconcentration Factor

CAS Chemical Abstracts Service

CLP Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures

DNEL Derived no-effect level

EC Identification code for each substance listed in EINECS

EC₅₀ Concentration of a substance when it is affected 50% of the population

EINECS European Inventory of Existing Commercial Chemical Substances

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016	Version	2.0
Revision date	21. August 2019		

EmS	Emergency plan
EU	European Union
IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
IC ₅₀	Concentration causing 50% blockade
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry
LC ₅₀	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD ₅₀	Lethal dose of a substance in which it can be expected death of 50% of the population
LOAEC	Lowest observed adverse effect concentration
LOAEL	Lowest observed adverse effect level
log K _{ow}	Octanol-water partition coefficient
MARPOL	International Convention for the Prevention of Pollution From Ships
NOAEC	No observed adverse effect concentration
NOAEL	No observed adverse effect level
NOEC	No observed effect concentration
NOEL	No observed effect level
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, Bioaccumulative and Toxic
PEL	Permissible Exposure Limit
PNEC	Predicted no-effect concentration
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Agreement on the transport of dangerous goods by rail
UN	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very Persistent and very Bioaccumulative

Eye Dam. Serious eye damage

Skin Irrit. Skin irritation

Training guidelines

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

Recommended restrictions of use

Designed for professional use only.

Information about data sources used to compile the Safety Data Sheet

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. The Act No. 350/2011 Coll., on Chemical Substances and Chemical Preparations as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

The changes (which information has been added, deleted or modified)

2.0 General revision of all sections of the safety data sheet according to Regulation (EC) No 830/2016 of the European Parliament and of the Council

Statement

ERBA LYTE PRO PLUS REAGENT PACK_SOL R

Creation date	18. January 2016		
Revision date	21. August 2019	Version	2.0

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.