

Product Guide

HOMOCYSTEINE 2R 200

Pack name	FHRWR 200	MSDS
Cat. No.	50003782	
Reagent 1	1 x 60 ml	MSDS HOMOCYSTEINE 2R 200_R1
Reagent 2	1 x 10 ml	<i>Note 1:</i> MSDS HOMOCYSTEINE 2R 200_R2
Calibrator 1 (CAL)	1 x 3 ml	<i>Note 2</i>
Calibrator 2 (CAL)	1 x 3 ml	<i>Note 2</i>

Calibrator 1

Homocystein blank (0 µmol/l)

Calibrator 2

Homocystein solution (28 µmol/l)

Note 1: Safety data sheet is compiled at a request of the recipient (in compliance of Article 31(3) of European Regulation (EC) No 1907/2006 (REACH)).

Note 2: The reagent is not classified as dangerous according to Regulation (EC) No 1272/2008 (CLP). In compliance with European Regulation (EC) No 1907/2006 (REACH) and harmonized Czech laws, safety data sheet for this product is not necessary.

Name of manufacturer:

Erba Lachema s.r.o.

Place of business:

Brno, Karásek 2219/1d

Postcode 621 00

Czech Republic

Phone:

+420 517 077 111

E-mail:

msds@erba.com