

ADVISORY NOTICE

02_2022_EC90 SW version 1.10.1 rev. 2


Dear customer,

We hereby inform you, that the device software of EC90 must be updated to SW version 1.10.1 (R1.01-V0.22) issued 20 January 2022. The SW version 1.10.1 is already included from device produced from 24th January 2022 onwards


For SW version 1.10.1, please contact your local service engineer or distributor to set up new SW version from web side [ERBA Academy](#).

To highlight few commonly asked questions we are adding few information points:

- Only FAT/FAT32 USB drives are compatible with EC 90 instrument
- Printout is customizable. These are options:



- Default setting for QC strictness is **None**, however we can only recommend to operate the instrument with higher QC strictness.



Applied for:
INSTRUMENT:

EC 90 Analyser

cat.no.INS00079

Erba Lachema s.r.o.



In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management

A handwritten signature in black ink, consisting of a stylized first name and a last name.

Graz, 21.01.2022

Brno 7th February 2022

Ref.No.AN_02_2022

