

ADVISORY NOTICE

URIC ACID 275 Cat. No.: XSYS0021

Dear customers,

We would like to draw your attention to the correction of the application for the UA 275 diagnostic kit on XL automatic analyzers:

- Change only applies to application using SI units ($\mu\text{mol/l}$)
- Corrected Assay type: 2-point (instead of former 1-point)

Please, fix this application in your XL analyzer.

In case of problems or questions, please do not hesitate to contact Erba Lachema Service support.

Best regards,

Erba Lachema s.r.o.
Technical product management

Brno 10th June 2022

