

ADVISORY NOTICE

Change of EC cartridge packaging

Dear valued customer,

To help provide better environmental control of EC Cartridge Biosensors during transportation, we would like to inform you that from 27th June 2022, all biosensors will be individually foil packed separately from the cartridge.

To prepare the cartridge for use, all that is required is to remove the EC Biosensor from the foil and insert into EC cartridge as described in the user manual.

We thank you in advance for your understanding – if you have any questions, please contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.

Technical Product Management

Brno, 23rd July 2022

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

