

ADVISORY NOTICE

23_2022_EC90 SW version 1.12.3

Dear customer,

We hereby inform you, that the device software of EC90 is highly recommended for update to SW version 1.12.3 issued on 14th October 2022. The SW version 1.12.3 is already included in devices produced from SN 4664-2-202209 onwards.

The software changes include the fixing of known bugs as well as implementing of new features (among others):

- Improved RFID detection method
- High precision mode is removed
- Bio Sensor Lot Number is now 6 or 7 digits long
- Conductivity limit is unified for Cal1 and Cal2 solutions
- Improvement of automatic calibration cycle in case of Adjustment Curve Error
- Ukrainian Keyboard is now fully supported

For SW version 1.12.3, please contact your local service engineer or distributor to set up new SW version from web side [ERBA Academy](#).

Applied for:

INSTRUMENT:

EC 90 Analyser

cat.no.INS00079

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management



Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-

Brno 19th October 2022

Ref.No.AN_23_2022



Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, Czech Republic
Identification number: 269 18 846. Tax identification number: CZ26918846
Incorporated in the Commercial Register maintained by the Regional Court in Brno, Section C, insert 45458
Tel.: +420 517 077 111, e-mail: diagnostics@erbamannheim.com, www.eralachema.com