

## URINORM XL

Creation date	26th February 2021	Version	2.0
Revision date	21st May 2024		

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Substance / mixture	URINORM XL
Number	mixture
Other mixture names	REG00060
	URINORM N XL; URINORM P XL

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Mixture's intended use**

URINORM XL control urines are liquid synthetic material with concentrations in the normal (URINORM XL N) and pathological (URINORM XL P) ranges. URINORM XL control urines are designed to verify the precision and accuracy of PHAN AUTO diagnostic strips and automatic urine analyzer LAURA XL intended for in vitro diagnostics. This product is intended also for evaluation of sediment. For professional use only.

**Main intended use**

PC-MED-OTH Other medical devices

**Secondary uses**

PC-TEC-19 Reagents and laboratory chemicals

**Mixture uses advised against**

not available

**1.3. Details of the supplier of the safety data sheet****Manufacturer**

Name or trade name	Erba Lachema s.r.o.
Address	Karásek 2219/1d , Brno, 62100
	Czech Republic
Identification number (CRN)	26918846
VAT Reg No	CZ26918846
Phone	+420 517 077 111
E-mail	msds@erba.com
Web address	www.erbalachema.com

**Competent person responsible for the safety data sheet**

Name	Erba Lachema s.r.o.
E-mail	msds@erba.com

**1.4. Emergency telephone number**

European emergency number: 112

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Classification of the mixture in accordance with Regulation (EC) No 1272/2008**

The mixture is not classified as dangerous according to Regulation (EC) No 1272/2008.

**2.2. Label elements**

none

**2.3. Other hazards**

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605. Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended.

## URINORM XL

Creation date 26th February 2021  
Revision date 21st May 2024 Version 2.0

## SECTION 3: Composition/information on ingredients

## 3.2. Mixtures

## Chemical characterization

Mixture of substances and additives specified below.

**Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment**

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note
CAS: 7647-14-5 EC: 231-598-3	Sodium chloride	<3		
Index: 017-002-01-X CAS: 7647-01-0 EC: 231-595-7	Hydrochloric acid 35%	<1	Met. Corr. 1, H290 Skin Corr. 1B, H314 STOT SE 3, H335 Specific concentration limit: Skin Corr. 1B, H314: $C \geq 25\%$ Skin Irrit. 2, H315; Eye Irrit. 2, H319: $10\% \leq C < 25\%$ STOT SE 3, H335: $C \geq 10\%$	1
Index: 011-002-00-6 CAS: 1310-73-2 EC: 215-185-5	sodium hydroxide	<0,1	Skin Corr. 1A, H314 Specific concentration limit: Skin Corr. 1B, H314: $2\% \leq C < 5\%$ Skin Corr. 1A, H314: $C \geq 5\%$ Eye Irrit. 2, H319: $0.5\% \leq C < 2\%$ Skin Irrit. 2, H315: $0.5\% \leq C < 2\%$	1
CAS: 57-88-5 EC: 200-353-2	Cholesterol	<0,1	not classified as dangerous	
Index: 011-004-00-7 CAS: 26628-22-8 EC: 247-852-1	sodium azide	$\leq 0,09$	Acute Tox. 2, H300 Aquatic Acute 1, H400 (M=1) Aquatic Chronic 1, H410 (M=1) EUH032	1
Index: 603-117-00-0 CAS: 67-63-0 EC: 200-661-7	isopropanol	<0,01	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336	1

## Notes

1 A substance for which exposure limits are set.

Full text of all classifications and hazard statements is given in the section 16.

## SECTION 4: First aid measures

## 4.1. Description of first aid measures

Take care of your own safety. If any health problems are manifested or if in doubt, inform a doctor and show him information from this safety data sheet.

## If inhaled

Terminate the exposure immediately; move the affected person to fresh air. Seek medical attention if breathing difficulty occurs.

## If on skin

Remove contaminated clothes. Seek medical attention if symptoms of skin irritation appear.

## If in eyes

Rinse eyes immediately with a flow of running water, open the eyelids (also using force if needed); remove contact lenses immediately if worn by the affected person.

## If swallowed

Rinse out the mouth with clean water. In the event of issues, find medical help.

## URINORM XL

Creation date	26th February 2021	Version	2.0
Revision date	21st May 2024		

**4.2. Most important symptoms and effects, both acute and delayed****If inhaled**

Not expected.

**If on skin**

Not expected.

**If in eyes**

Not expected.

**If swallowed**

Not expected.

**4.3. Indication of any immediate medical attention and special treatment needed**

Symptomatic treatment.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**

Accommodate extinguishing components to the location of fire.

**Unsuitable extinguishing media**

not available

**5.2. Special hazards arising from the substance or mixture**

Inhalation of hazardous degradation (pyrolysis) products may cause serious health damage.

**5.3. Advice for firefighters**

Self-Contained Breathing Apparatus (SCBA) with chemical resistant gloves. Use a self-contained breathing apparatus and full-body protective clothing.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Follow the instructions in the Sections 7 and 8. Ensure adequate ventilation. Handle as potentially infectious.

**6.2. Environmental precautions**

Prevent contamination of the soil and entering surface or ground water.

**6.3. Methods and material for containment and cleaning up**

Wipe up with adsorbent material. Place in suitable container for prompt disposal. Label the container as to the potential infectious hazard. Spill areas can be decontaminated with 0.5% sodium hypochlorite, e.g., a fresh 1:10 dilution of common household bleach.

**6.4. Reference to other sections**

See the Section 7, 8 and 13.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling**

Prevent formation of gases and vapours in concentrations exceeding the occupational exposure limits. Use personal protective equipment as per Section 8. Observe valid legal regulations on safety and health protection. Wear protective safety glasses, gloves and clothing.

Wear disposable gloves while handling reagent. Thoroughly wash hands afterwards.

**7.2. Conditions for safe storage, including any incompatibilities**

Store in tightly closed containers in cold, dry and well ventilated areas designated for this purpose.

Storage temperature

min 2 °C, max 8 °C

**7.3. Specific end use(s)**

For in vitro diagnostic devices.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters**

The mixture contains substances for which occupational exposure limits are set.

## URINORM XL

Creation date

26th February 2021

Revision date

21st May 2024

Version

2.0

## Czech Republic

## Government Regulation 330/2023 Coll.

Substance name (component)	Type	Value	Note
Hydrochloric acid 35% (CAS: 7647-01-0)	PEL	8 mg/m <sup>3</sup>	irritating to mucous membranes (eyes, respiratory system) and skin
	PEL	5 ppm	
	NPK-P	15 mg/m <sup>3</sup>	
	NPK-P	10 ppm	
sodium hydroxide (CAS: 1310-73-2)	PEL	1 mg/m <sup>3</sup>	irritating to mucous membranes (eyes, respiratory system) and skin
	NPK-P	2 mg/m <sup>3</sup>	
sodium azide (CAS: 26628-22-8)	PEL	0,1 mg/m <sup>3</sup>	skin penetration is significantly involved during exposure, irritating to mucous membranes (eyes, respiratory system) and skin
	NPK-P	0,3 mg/m <sup>3</sup>	
isopropanol (CAS: 67-63-0)	PEL	500 mg/m <sup>3</sup>	irritating to mucous membranes (eyes, respiratory system) and skin
	PEL	200 ppm	
	NPK-P	1000 mg/m <sup>3</sup>	
	NPK-P	400 ppm	

## European Union

## Commission Directive 2000/39/EC

Substance name (component)	Type	Value	Note
Hydrochloric acid 35% (CAS: 7647-01-0)	OEL 8 hours	8 mg/m <sup>3</sup>	
	OEL 8 hours	5 ppm	
	OEL 15 minutes	15 mg/m <sup>3</sup>	
	OEL 15 minutes	10 ppm	
sodium azide (CAS: 26628-22-8)	OEL 8 hours	0,1 mg/m <sup>3</sup>	Skin
	OEL 15 minutes	0,3 mg/m <sup>3</sup>	

## URINORM XL

Creation date	26th February 2021		
Revision date	21st May 2024	Version	2.0

**8.2. Exposure controls**

Do not eat, drink and smoke during work. Wash your hands thoroughly with water and soap after work and before breaks for a meal and rest.

**Eye/face protection**

Wear Safety Glasses

**Skin protection**

Wear disposable gloves. Thoroughly wash hands.

**Respiratory protection**

Not required with adequate ventilation.

**Thermal hazard**

Not available.

**Environmental exposure controls**

Observe usual measures for protection of the environment, see Section 6.2.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state	liquid
Colour	Straw colored
Odour	Mild odor
Melting point/freezing point	data not available
Boiling point or initial boiling point and boiling range	data not available
Flammability	The product is non-flammable.
Lower and upper explosion limit	data not available
Flash point	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
pH	data not available
Kinematic viscosity	data not available
Solubility in water	data not available
Partition coefficient n-octanol/water (log value)	data not available
Vapour pressure	data not available
Density and/or relative density	data not available
Relative vapour density	data not available
Particle characteristics	data not available

**9.2. Other information**

not available

**SECTION 10: Stability and reactivity****10.1. Reactivity**

not available

**10.2. Chemical stability**

The product is stable under normal conditions.

**10.3. Possibility of hazardous reactions**

Unknown.

**10.4. Conditions to avoid**

The product is stable and no degradation occurs under normal use. Protect against flames, sparks, overheating and against frost.

**10.5. Incompatible materials**

Protect against strong acids, bases and oxidizing agents.

**10.6. Hazardous decomposition products**

Not developed under normal uses.

## URINORM XL

Creation date 26th February 2021  
Revision date 21st May 2024 Version 2.0

## SECTION 11: Toxicological information

## 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

No toxicological data is available for the mixture.

**Acute toxicity**

Based on available data the classification criteria are not met.

sodium azide					
Route of exposure	Parameter	Value	Exposure time	Species	Sex
Oral	LD <sub>50</sub>	27 mg/kg bw			
Inhalation	LC <sub>50</sub>	54 mg/m <sup>3</sup>	4 hours	Rat	

**Skin corrosion/irritation**

Based on available data the classification criteria are not met.

**Serious eye damage/irritation**

Based on available data the classification criteria are not met.

**Respiratory or skin sensitisation**

Based on available data the classification criteria are not met.

**Germ cell mutagenicity**

Based on available data the classification criteria are not met.

**Carcinogenicity**

Based on available data the classification criteria are not met.

**Reproductive toxicity**

Based on available data the classification criteria are not met.

**Toxicity for specific target organ - single exposure**

Based on available data the classification criteria are not met.

**Toxicity for specific target organ - repeated exposure**

Based on available data the classification criteria are not met.

**Aspiration hazard**

Based on available data the classification criteria are not met.

## 11.2. Information on other hazards

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

## SECTION 12: Ecological information

## 12.1. Toxicity

not available

**Acute toxicity**

Hydrochloric acid 35%				
Parameter	Value	Exposure time	Species	Environment
LC <sub>50</sub>	20.5 mg/l	96 hours	Fish	

## URINORM XL

Creation date 26th February 2021  
Revision date 21st May 2024 Version 2.0

sodium azide				
Parameter	Value	Exposure time	Species	Environment
LC <sub>50</sub>	680 µg/l		Fish	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	400 µg/l		Invertebrates	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	150 µg/l		Invertebrates	Salt water
EC <sub>50</sub> /LC <sub>50</sub>	348 µg/l		Algae	Fresh water
EC <sub>50</sub> /LC <sub>50</sub>	5.6 mg/l		Microorganisms	
NOEC	30 µg/l		Microorganisms	

**12.2. Persistence and degradability**

not available

**12.3. Bioaccumulative potential**

Not available.

**12.4. Mobility in soil**

Not available.

**12.5. Results of PBT and vPvB assessment**

Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

**12.6. Endocrine disrupting properties**

The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605.

**12.7. Other adverse effects**

not available

**SECTION 13: Disposal considerations****13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Proceed in accordance with valid regulations on waste disposal. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Empty containers may be used at waste incinerators to produce energy or deposited in a dump with appropriate classification. Perfectly cleaned containers can be submitted for recycling.

**Waste management legislation**

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

**SECTION 14: Transport information****14.1. UN number or ID number**

not subject to transport regulations

**14.2. UN proper shipping name**

not relevant

**14.3. Transport hazard class(es)**

not relevant

**14.4. Packing group**

not relevant

**14.5. Environmental hazards**

not relevant

**14.6. Special precautions for user**

Reference in the Sections 4 to 8.

**14.7. Maritime transport in bulk according to IMO instruments**

not relevant

## URINORM XL

Creation date	26th February 2021		
Revision date	21st May 2024	Version	2.0

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. Commission Regulation (EU) 2020/878 of 18 June 2020 amending Annex II to Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

**National regulations (Germany)**

WGK Water hazard class: WGK 1 - slightly hazardous to water

**15.2. Chemical safety assessment**

not available

**SECTION 16: Other information****A list of standard risk phrases used in the safety data sheet**

H225	Highly flammable liquid and vapour.
H290	May be corrosive to metals.
H300	Fatal if swallowed.
H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.

**A list of additional standard phrases used in the safety data sheet**

EUH032	Contact with acids liberates very toxic gas.
--------	--

**Other important information about human health protection**

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

**Key to abbreviations and acronyms used in the safety data sheet**

ADR	European agreement concerning the international carriage of dangerous goods by road
BCF	Bioconcentration Factor
CAS	Chemical Abstracts Service
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures
EC	Identification code for each substance listed in EINECS
EINECS	European Inventory of Existing Commercial Chemical Substances
EmS	Emergency plan
EU	European Union
EuPCS	European Product Categorisation System
IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry



## URINORM XL

Creation date	26th February 2021	Version	2.0
Revision date	21st May 2024		

LC <sub>50</sub>	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD <sub>50</sub>	Lethal dose of a substance in which it can be expected death of 50% of the population
log K <sub>ow</sub>	Octanol-water partition coefficient
NOEC	No observed effect concentration
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, Bioaccumulative and Toxic
PEL	Permissible Exposure Limit
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Agreement on the transport of dangerous goods by rail
UN	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very Persistent and very Bioaccumulative
Acute Tox.	Acute toxicity
Aquatic Acute	Hazardous to the aquatic environment
Aquatic Chronic	Hazardous to the aquatic environment (chronic)
Flam. Liq.	Flammable liquid
Met. Corr.	Corrosive to metals
Skin Corr.	Skin corrosion
STOT SE	Specific target organ toxicity - single exposure

**Training guidelines**

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

**Recommended restrictions of use**

not available

**Information about data sources used to compile the Safety Data Sheet**

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended.  
REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

**The changes (which information has been added, deleted or modified)**

2.0 General revision of all sections of the safety data sheet according to Commission Regulation (EU) 2020/878.

**More information**

Classification procedure - calculation method.

**Statement**

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.