

AN 2022\_34

## ADVISORY NOTICE

### XL MULTICAL

Dear customers,

With a mission to continuously enhance user experience and satisfaction, we would like to inform you regarding the packaging change and the related pack name change of biochemistry calibrator XL MULTICAL; a larger pack containing 10 vials will also be available – please see table below:

CURRENT STATUS:			NEW PACKAGING:			
Cat. No.:	Name/Pack name:	Content:	Cat. No.:	Name:	Pack name:	Content:
XSYS0034	XL MULTICAL	4x3 mL	XSYS0034	XL MULTICAL	XL MULTICAL 4x3	4x3 mL
			XSYS0122	XL MULTICAL	XL MULTICAL 10x3	10x3 mL

In connection with this change, the design of value sheet was modified and several values for new analytes were added:

- Alkaline phosphatase – DGKC method with MEG buffer
- Alanine aminotransferase – IFCC method with PDP
- Aspartate aminotransferase – IFCC method with PDP
- Amylase total – IFCC method
- Amylase pancreatic – method EPS with maltoheptaoside
- LDH – IFCC method (L to P)
- Uric acid – colorimetric method with AOD

Please note, that upcoming change has no influence on chemical composition, **the matrix of biochemistry calibrator remains unchanged.**

In case of any questions, please do not hesitate to contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.  
Technical Product Management

Brno, 19<sup>th</sup> December 2022

