

ADVISORY NOTICE

36_2022_LAURA XL SW version 1.2.1

Dear customer,

We hereby inform you, that the device software of LAURA XL is highly recommended for update to SW version 1.2.1 issued 20th December 2022. The SW version 1.2.1 is already included from device produced from SN 0364-X2-180822 onwards.

The new SW version 1.2.1 includes a correction of a bug which caused the non-transmission by LIS connection of chemistry results when the strips used are UNDEKAPHAN auto.

User manual for SW version 1.2.0 is also valid for the new SW version 1.2.1.

For SW version 1.2.1, please contact your local service engineer or distributor to set up new SW version from web side [ERBA Academy](#).

Applied for:

INSTRUMENT:



LAURA XL

cat.no.INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -22-

Brno 20th December 2022

Ref. No.AN_36_2022