

AN 2023_05

ADVISORY NOTICE

Update of target ranges for Bio-Rad Liquichek – batches: 87740 and 87750

LAURA XL, Catalogue No.: INS00065
Bio-Rad Liquichek Urinalysis Control, Catalogue No.: 435

Dear customers,

we would like to inform you, that according to our latest observation and customers feedback we are forced to update target ranges of UBG and PRO zones valid for LAURA XL instrument combined with DEKAPHAN AUTO diagnostic strips for Bio-Rad Liquichek Urinalysis Control – batches 87740 (exp: 07/2023), 87750 (exp: 12/2023).

According to our additional testing, the updated values are following:


Bio-Rad Liquichek Urinalysis Control, batch 87740 (exp: 07/2023)		
Analyte	Level 1 – 87741	Level 2 – 87742
Protein	negative	30–500 mg/dl (0.3–5 g/l) (1+ - 3+)
Urobilinogen	normal	3–12 mg/dl (51–203 µmol/l) (2+ - 4+)

Bio-Rad Liquichek Urinalysis Control, batch 87750 (exp: 12/2023)		
Analyte	Level 1 – 87751	Level 2 – 87752
Protein	negative	30–500 mg/dl (0.3–5 g/l) (1+ - 3+)
Urobilinogen	normal	3–12 mg/dl (51–203 µmol/l) (2+ - 4+)

In case of any questions, please contact your local ERBA team.

Best regards,

Erba Lachema s.r.o.
Technical Product Management
Brno, 14th February 2023


Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-