

AN_2023_16

ADVISORY NOTICE

MIKROLATEST Beta-LACTAMtest
Catalog No. MLT00034

Dear customer,

We would like to inform you about a change in the instructions for the β -LACTAM test. Based on EUCAST recommendations, the test is no longer suitable for the detection of *Haemophilus* as stated in the original instructions.

The β -LACTAM test is intended for the rapid acidimetric detection of bacterial β -lactamase activity. This acidimetric method is recommended for use in testing β -lactamase production by *Staphylococcus* species and *Neisseria*.

If you have any questions, please do not hesitate to contact us.

Best regards

Erba Lachema s.r.o.
Technical Product Management
Brno, 21.06.2023

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

