

AN 2023_21

ADVISORY NOTICE

21_2023_H360 cat.no INS00077

Dear customer,

We hereby inform you, that for device H360 cat.no INS00077 new update was issued in January 2023 for instruments H360 (since SNs K1001224xxxx).
Mentioned update include HW changes and SW changes which are related to the HW changes.

- 1) HW changes
 - 1.a Main Board_ 22.01.0288A (separate the FPGA chip)
 - 1.b Pressure sensor_ 22.01.0022B third generation
- 2) The new SW version for device H360 is **A14.10**

This new change will not affect the clinical performance, customer usage and it is compatible with the old type analyser. Hence, there will be no risk involved.

The SW version A14.10 was placed on our web side, please visit [ERBA Academy](#).

In case you have any additional question, please do not hesitate to contact us.

Best regards

Erba Lachema s.r.o.
Technical Product Management

Brno 24th July 2023



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Ref.No.AN_21_2023