

AN 2023_22

ADVISORY NOTICE

22_2023_CHEM-7 cat.no INS00014

Dear customer,

We hereby inform you, that for device CHEM-7 cat.no INS00014 new update was issued in January 2023 for instruments CHEM-7 (since SNs S30xx).
Mentioned update include HW changes and SW changes which are related to the HW changes.

1) HW changes

- ASSEMBLY PCB Arm Controller CPU EL20205R3.1 with microcontroller LPC1788 This Assembly PCB Arm Controller CPU EL20205R3.1 is compatible with board M3 PCB EL20205R3.1 with NEW ADC PCB & OLD ADC PCB (ELLC2002R2) with CP2102 drivers to be installed for

2) The new SW version for device CHEM-7 is SW 2022.C0.09

- Bug fixed and improvements
- Embedded SW for mainboard cortex EL20205R3.1 (Board with Microcontroller LPC1788)

This new change will not affect the clinical performance, customer usage.
Hence, there will be no risk involved.

The SW version SW 2022.C0.09 was placed on our web side, please visit [ERBA Academy](#).

In case you have any additional question, please do not hesitate to contact us.

Best regards

Erba Lachema s.r.o.
Technical Product Management



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno 24th July 2023

Ref.No.AN_22_2023