

AN 2023_24

ADVISORY NOTICE

24_2023_XL180 cat.no INS00037 and INS00038

Dear customer,

We hereby inform you, that for device XL180 cat.no INS00037 and INS00038 new update was issued in February 2023, instruments XL180_SN 160681 onwards.

Application software ver.XL180 SW V2022.C02_V2022.R.01 is released to include new feature, improvements & to address issues reports in previous software version.

The new SW versions for device **XL180 cat.no INS00037 and INS00038_ASW XL180 _V2022.C02** along with embedded software **ESW ver. XL180_V2022.R.01** include changes:

- 1) Master Calibration - Provision to load Master Calibration from RF-Card
- 2) Provision for Quality Control of Urine Parameter, Unit Conversion
- 3) PCR TD/RS/010: Rgt Abs Min/Max check for Flagging Calibration.
- 4) PCR PG/RS/007: Absorbance above 3.0 – Flag with Result NA and facility to Re-run.
- 5) Facility to user to edit Decimal Places in ISE Test Parameters (similar to Photometric Tests).
- 6) Database Security Enhanced.
- 7) Reagent Bottle Request/Transfer facility in Consumable Management Utility
- 8) EL's AS-42 (issue 13-Sep-2021, Conventional unit) pre-loaded in Blank Database..

This new change will not affect the clinical performance, customer usage.
Hence, there will be no risk involved.

The SW versions V2022.C02_V2022.R.01 were placed on our web side, please visit [ERBA Academy](#).

In case you have any additional question, please do not hesitate to contact us.

Best regards

Erba Lachema s.r.o.
Technical Product Management



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno 24th July 2023

Ref.No.AN_24_2023

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, Czech Republic
Identification number: 269 18 846. Tax identification number: CZ26918846
Incorporated in the Commercial Register maintained by the Regional Court in Brno, Section C, insert 45458
Tel.: +420 517 077 111, e-mail: diagnostics@erba.com, www.erbalachema.com