

**ADVISORY NOTICE**

27\_2023\_CRP CONTROL values

Product name: CRP CONTROL HIGH, Cat. No. BLT20013  
CRP CONTROL LOW, Cat. No. BLT20014

Dear customer,

We would like to inform you that the value sheet for the mentioned products has been extended with values for XL 200, XL 640, and XL 1000.

Since the new production lot, the values will be stated in the Instruction for Use attached to the kit.

The values for the lots already on the market are as follows:

CRP CONTROL HIGH: lot no 2301086 and lot no 2305151

Analyser	Value	Interval	1 SD	Unit
XL-200	3.74	2.99-4.49	0.25	mg/dl
	37.4	29.9-44.9	2.49	mg/l
XL-640	3.61	2.89-4.33	0.24	mg/dl
	36.1	28.9-43.3	2.41	mg/l
XL-1000	3.34	2.68-4.00	0.22	mg/dl
	33.4	26.8-40.0	2.23	mg/l

CRP CONTROL LOW: lot no 2301103 and lot no 2305152

Analyser	Value	Interval	1 SD	Unit
XL-200	1.63	1.30-1.96	0.11	mg/dl
	16.3	13.0-19.6	1.09	mg/l
XL-640	1.63	1.30-1.96	0.11	mg/dl
	16.3	13.0-19.6	1.09	mg/l
XL-1000	1.56	1.26-1.86	0.10	mg/dl
	15.6	12.6-18.6	1.04	mg/l

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.  
Technical product management

Brno 22<sup>nd</sup> August 2023

**Erba Lachema s.r.o.**  
Karásek 2219/1d, 621 00 Brno  
DIČ: CZ26918846 -7-

Ref.No.AN 2023\_27