

## ADVISORY NOTICE

26\_2023\_LAURA XL SW version 1.2.3.31

Dear customer,

We hereby inform you that the updated software (version 1.2.3.31) for LAURA XL automatic urine analyser was released on the 14<sup>th</sup> of August 2023. The updated SW version includes mainly repair of minor bugs detected by manufacturer including customers' requests from global market, compatibility with new 16-way valve, improved hydraulic system and addition of UA language into user interface.

Erba Lachema highly recommend you update the SW on your instrument to the latest version (1.2.3.31).

For the latest SW version, please contact your local service engineer or distributor to install updated SW version to your instrument. The updated SW version you can find on ERBA Academy website.

Applied for:

**INSTRUMENT:**

LAURA XL

cat.no.: INS00065

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.  
Technical product management



**Erba Lachema s.r.o.**  
Karásek 2219/1d, 621 00 Brno  
DIČ: CZ26918846 -7-

Brno 24<sup>th</sup> August 2023

Ref. No.AN\_26\_2022