

ADVISORY NOTICE

30_2023_EC90 SW version 1.20.8

Dear customer,

We hereby inform you, that the device software of EC90 is mandatory for update to SW version 1.20.8 issued 22nd September 2023 due to RFID communication improvement. EC 90 analyser SW update from V1.13.3 to V1.20.8 implements new features based on the additional requirements featuring (among others):

- Manual triggered LIS data export; removing all GMS
- Measurement aspirating and cartridge emptying warning, pop up for the user
- Improved cleaning routine
- Biosensor on board stability exceeded with grey parameter button on the main screen
- RFID communication Improved resolving 0 sample 0 days challenge with valid cartridges

The SW version 1.20.8 is already included in devices produced from SN 8064-2-202308 onwards.

For SW version 1.20.8, please contact your local service engineer or distributor to set up new SW version from web side [ERBA Knowledge Base](#).

Applied for:

INSTRUMENT:

EC 90 Analyser

cat.no.INS00079

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management



Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-

Brno 3th October 2023

Ref.No.AN_30_2023

