

ADVISORY NOTICE

33_2023_ change of color of negative control HIV Ab&Ag kit

Product name: HIV Ab&Ag kit

Cat. No. IVCOMB.CE

Dear customer,

We would like to inform you that the color of the Negative Control of the HIV Ab&Ag kit (code IVCOMB.CE) and of any other Elisa kit produced by Diapro can vary from yellowish to brownish depending on the LOT of serum or defibrinated plasma used.

The performance of the Negative Control always complies with the assigned requirements (in the case of HIV Ab&Ag OD < 0.200 after blanking), regardless of the color tonality.

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Nemcova!

Erba Lachema s.r.o.
Technical product management

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-

Brno 19th October 2023

Ref.No.AN 2023_33