

ADVISORY NOTICE

37_2023_Shock sensors on Instrument packaging


Dear customer,

We hereby inform you, that some Erba Lachema instruments in transport box are shipped with Shock sensors. The function of the sensors is to inform the end user, whether the manipulation with instrument was according to specified transport conditions.

If the sensors are activated, that means there might be an impact to instrument functionality. In case transport box contains Shock sensor, please check the sensor within a delivery of instrument. Once instrument is accepted on the stock it is a confirmation there is no packaging damage or shock sensor problem.

If the sensors are activated, keep evidence and contact delivery company.

Example of Shock sensors:



In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management

Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno 09th November 2023

Ref.No.AN_37_2023

