

ADVISORY NOTICE

07_2024_ Erba Thrombin Reagent

Dear customer,

We hereby inform you that if you are using Erba Thrombin Reagent 10/20 on Erba ECL semi-automated analysers ECL105/412 or automated analysers ECL 760 together with control material Erba Control N (lot **2401077/2401088**), following values assigned to this reagent-control lot combination should be used:

Product/Test Producto/Prueba	Unit Jednotka Unidad	ECL 105 / 412			ECL 760		
		Optical Optická detekce Detección óptica			Optical Optická detekce Detección óptica		
		Value Hodnota Valor	Range Interval Rango	SD value	Value Hodnota Valor	Range Interval Rango	SD value
FIBRINOGEN / Erba Thrombin Reagent (for determination of Fibrinogen) / Clauss Fibrinogen	g/l	3.13	2.66–3.59	0.16	3.61	3.07–4.15	0.17

Following Erba products are concerned:

Erba Thrombin Reagent 10		Cat. no. EHL00005
Erba Thrombin Reagent 20		Cat. no. EHL00006
Erba Control N	Lot 2401077/2401088	Cat. no. EHL00014

When used on:

ECL105 semi-automated coagulometer	Cat. no. INS00059
ECL412 semi-automated coagulometer	Cat. no. INS00060
ECL760 automated coagulometer	Cat. no. INS00070

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management**Erba Lachema s.r.o.**
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 2.5.2024

Ref. No.AN_07_2024