

ADVISORY NOTICE

10_2024_Medica Corporation – electrodes (shelf life/onboard lifetime)

Dear customer,

We would like to provide you with clear explanation of the “expiration” or “install by” date and “onboard lifetime” of the electrodes produced by Medica Corporation.

The “expiration” or “install by” date on the package indicates the product's shelf life. The onboard lifetime begins on the installation date, if the electrode is installed at any time within ‘shelf life’ period.

Examples:

- 1.) If the expiration/install-by date is May 1st, 2024, but the product was installed in September 2023, the 6 months period started in September 2023 already.
- 2.) If the expiration/install-by date is May 1st, 2024, and the product is installed on April 30th, 2024, the 6 months period starts on the installation date (in this case April 30th, 2024).

Affected products

Catalogue number	Part Number	Product name	Shelf Life	Onboard lifetime
5201	51000144	ISE NA ELECTRODE FOR ERBA XL	1 year	6 months/10000 samples
5202	51000145	ISE K ELECTRODE FOR ERBA XL	1 year	6 months/10000 samples
5204	51001521	ISE REFERENCE ELECTRODE FOR ERBA XL 113952	1 year	6 months/10000 samples
5205	51001519	ISE LI ELECTRODE FOR ERBA XL	1 year	6 months/10000 samples
5207	51000143	ISE CL ELECTRODE FOR ERBA XL	1 year	6 months/10000 samples

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.
Technical product management



Erba Lachema s.r.o.
Karásek 2219/1d, 621 00 Brno
DIČ: CZ26918846 -7-

Brno, 17th May 2024

Ref. No.AN_10_2024

